
**THUASNE**

# Patella Reliever

<b>fr</b>	Genouillère articulée avec système de réalignement de la rotule	6
<b>en</b>	Hinged knee orthosis with patella realignment system	7
<b>de</b>	Gelenkorthese mit System zur Ausrichtung der Kniescheibe	9
<b>nl</b>	Kniebrace met scharnierspalken en patella-uitlijningsstelsel	10
<b>it</b>	Ginocchiera articolata con sistema di riallineamento della rotula	12
<b>es</b>	Rodillera articulada con sistema de realineamiento de la rótula	14
<b>pt</b>	Joelheira articulada com sistema de realinhamento da rótula	15
<b>da</b>	Artikuleret knæbind med system til at tilpasse knæskallens position	17
<b>fi</b>	Nivelletty polvituki patellan uudelleenkohdistusjärjestelmä	19
<b>sv</b>	Ledat knästöd med justeringssystem för knäskålen	20
<b>el</b>	Κνηδεμόνας γόνατος με ενσωματωμένες αρθρώσεις και σύστημα έλξης της επιγονατίδας	22
<b>cs</b>	Kloubová ortéza kolenního kloubu se systémem vyrovnání číšky	23
<b>pl</b>	Zginalna orteza stawu kolanowego z systemem wyrównania rzepki	25
<b>lv</b>	Artikulēta ceļa ortoze ar ceļa skriemeļa izlīdzināšanas sistēmu	26
<b>lt</b>	Lankstus antkelis su girtelės sulyginimo sistema	28
<b>et</b>	Liigendatud põlvetugi põlvekedra korrigeerimise süsteemiga	29
<b>sl</b>	Zgibna opornica s sistemom za poravnavo pogačice	31
<b>sk</b>	Kolenná ortéza s patelárnym trakčným systémom	32
<b>hu</b>	Ízületes térdortézis térdkalács-beállító rendszerrel	34
<b>bg</b>	Наколенка със система за изправяне на пателата	36
<b>ro</b>	Genunchieră cu articulație și sistem de realiniere a rotulei	37
<b>ru</b>	Шарнирная повязка для коленного сустава с системой настройки для коленной чашечки	39
<b>hr</b>	Ortoza za koljenski zglob sa sustavom za podešavanje položaja patele	41
<b>zh</b>	带有调整髌骨系统的关节护膝	42
<b>ar</b>	دعامة ركبة مفصليّة مع نظام محاذاة للرصفّة	45

	1	32 - 34 cm
	2	35 - 37 cm
	3	38 - 41 cm
	4	42 - 45 cm
	5	46 - 49 cm
	6	50 - 53 cm

<b>fr</b>	Circonférence du genou
<b>en</b>	Knee circumference
<b>de</b>	Knieumfang
<b>nl</b>	Omvang van de knie
<b>it</b>	Circonferenza del ginocchio
<b>es</b>	Circunferencia de la rodilla
<b>pt</b>	Circunferência do joelho
<b>da</b>	Omkreds af knæ
<b>fi</b>	Polven ympärys
<b>sv</b>	Omkrets runt knä
<b>el</b>	Περίμετρος του γόνατος
<b>cs</b>	Obvod kolena
<b>pl</b>	Obwód kolana
<b>lv</b>	Celgala apkārtmērs
<b>lt</b>	Kelio apimtis
<b>et</b>	Põlve ümbermõõt
<b>sl</b>	Obseg kolena
<b>sk</b>	Obvod kolena
<b>hu</b>	Térd körfogata
<b>bg</b>	Обиколка на коляното
<b>ro</b>	Circumferința genunchiului
<b>ru</b>	Диаметр окружности колена
<b>hr</b>	Opseg koljena
<b>zh</b>	膝关节周长
<b>ar</b>	محيط الركبة

Correction biomécanique	Compression	Effet proprioceptif
Biomechanical correction	Compression	Proprioceptive effect
Biomechanische Korrektur	Kompression	Propriozeptive Wirkung
Biomechanische correctie	Compressie	Proprioceptieve werking
Correzione biomeccanica	Compressione	Effetto proprioceetivo
Corrección biomecánica	Compresión	Efecto propioceptivo
Correção biomecânica	Compressão	Efeito proprioceetivo
Biomekanisk korrektion	Kompression	Proprioceptisk effekt
Biomekaaninen korjaus	Kompressio	Liikeaistiin kohdistuva vaikutus
Biomekanisk korrigerig	Kompression	Proprioceptiv effekt
Βιο-μηχανικός έλεγχος	Συμπίεση	Αποτέλεσμα ιδιοδεκτικότητας
Biomechanická korekce	Komprese	Proprioceptivní účinek
Korekcja biomechaniczna	Kompresja	Propriocepcja
Biomehāniskā korekcija	Kompresija	Proprioceptīvs efekts
Mechaninė biologinė korekcija	Kompresija	Propriocepcinis poveikis
Biomehaaniline korrektsioon	Kompressioon	Proprioitseptiivne toime
Popravek biomehanske motnje	Kompresija	Proprioceptivni učinek
Biomechanická korekcia	Kompresia	Proprioceptívny účinok
Biomechanikus korrekció	Szorítás	Proprioceptív hatás
Биомеханична корекция	Компрессия	Проприоцептивный эффект
Corecție biomecanică	Compresie	Efect proprioceptiv
Биомеханическая коррекция	Компрессия	Проприоцептивный эффект
Biomehānička korekcija	Kompresija	Učínak proprioceptivnog vježbanja
生物力学矫正	压紧	本体感

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---

## en HINGED KNEE ORTHOSIS WITH PATELLA REALIGNMENT SYSTEM

### Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Product available in left or right version.

Available in 6 sizes.

The device is composed of:

- a half-moon shape patella insert ❶,
- a traction system for correction and support of the patella glide path ❷,
- an anatomically-shaped knitting integrating a light and soft knit at the back of the knee ❸,
- 2 rigid hinged uprights ❹,
- 1 inelastic strap ❺.

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### Properties/Mode of action

Biomechanical correction: half-moon shape patella insert with traction system to realign the patella.

An anatomically-shaped knitting for even compression.

The brace is maintained on the leg thanks to the silicone-coated threads at the top of the knee brace.

## Indications

- Femoro-patellar osteoarthritis (mild to moderate).
- Anterior knee or patella pain.
- Patellar instability.
- Patella luxation and subluxation.
- Abnormal patella displacement (patella alta).
- Femoro-patellar pain syndrome and/or patella instability.
- Post-operative use (after Medial Patellofemoral Ligament reconstruction, lateral release).

## Contraindications

- Do not use the product if the diagnosis has not been confirmed.
- Do not use in the event of known allergy to any of the components.
- Do not apply the product in direct contact with broken skin or an open wound without an adequate dressing.
- History of venous or lymphatic disorders.
- Do not use in the event of major venous thromboembolic history without thromboprophylaxis.

## Precautions

- Verify the product's integrity before every use.
- Do not use the device if it is damaged.
- Choose the appropriate size to fit the patient, referring to the size chart.
- It is recommended that a healthcare professional supervises the first application.
- Strictly comply with your healthcare professional's prescription and recommendations for use.
- This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.
- For hygiene and performance reasons, do not re-use the product for another patient.
- It is recommended to adequately tighten the device to achieve support/immobilisation without restricting blood circulation.
- In the event of discomfort, significant hindrance, pain, variation in limb volume, abnormal sensations or change in colour of the extremities, remove the device and consult a healthcare professional.
- In the event of any modification in the product's performance, remove it and consult a healthcare professional.
- Do not wear the product in a medical imaging machine.
- Do not use the device in case of application of certain products on the skin (creams, ointments, oils, gels, patches...).

## Undesirable side-effect

- This device can cause skin reactions (redness, itching, burns, blisters...) or wounds of various degrees of severity.
- Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

## Instructions for use/Application

- It is recommended to wear the device directly on the skin, unless contraindicated.
- Before you put it on loosen the strap.
  - To avoid damage to the knit fabric, please make sure that the self-gripping pads are always fastened to the designated self-gripping surfaces when putting it on and taking it off.
- Have the patient sit on the edge of a chair, and have the patient bend his/her knee to a position of approximately 30-60° of flexion.
- Slip the knee support on the leg with your foot extended.
  - To make it easier to put on, use the pull tabs at the top of the knee brace ❶.
- Make sure that your knee cap is positioned in the centre of the patella insert ❷.
- Close the calf strap ❸.
- To adjust the pressure of the patella traction system as needed, detach the strip of the self-fastening area on the interior side of the patella ❹. Pull the strip upward until you reach the desired pressure ❺. Refasten the strip on the self-fastening surface.
- To remove the knee brace: loosen the calf strap and pull on the lower part to remove it from the bottom.

## Care

- Close the self-fastening tabs before washing. Machine washable at 30°C (delicate programme).
- If possible use a washing net. Do not use detergents, fabric softeners or aggressive products

(products containing chlorine). Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). If the device is exposed to seawater or chlorinated water, make sure to rinse it in clear water and dry it.

#### **Storage**

Store at room temperature, preferably in the original packaging.

#### **Disposal**

Dispose of in accordance with local regulations.

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©Thuasne - 2026801 (2021-04)