




ThuasneUSA / Lantz Medical

*Phone: 866-236-8889


*Fax: 877-406-4872

*Web: www.thuasneusa.com

Stat-A-Dyne™ Pro/Sup *Fitting Instructions*



- Telescope the device to the appropriate length of the patient's forearm. Adjust the length by loosening the two screws on the bottom of the device. The length is adjusted correctly when the elbow is positioned 90°, while allowing room for the fingers to move freely.
- Once the device's length has been adjusted, apply the device by first putting the forearm into the distal cuff and inserting the thumb through the thumb hole. *The thumb should fit comfortably through the thumb hole avoiding pressure on the superficial branch of the radial nerve. The 2 screws on top of the forearm cuff are to allow adjustment for the thumb and fingers to fit comfortably.



- After the forearm is in the device, the humeral cuff may be pulled over the biceps/triceps. It is ideal to adjust the height of the humeral cuff so there are four finger widths between the top of the cuff and the crease of the axilla. The 2 screws posterior to the triceps are to allow both elevation / depression adjustments.
- Tighten and trim any excess length from the humeral and the forearm straps.
- The treatment session is now ready to begin.

Thanks for allowing us to be a part of your recovery!

Contact your Lantz Medical Sales Consultant for additional information on all our products or at www.thuasneusa.com