
SLEEQ EXT

INSTRUCTIONS FOR USE / GEBRAUCHSANWEISUNG
 INSTRUCCIONES DE USO
 INSTRUCTIONS D'UTILISATION
 ISTRUZIONI PER L'USO

1.0

en Start


es Empezar

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fr Démarrage

it Avvio

1.1

Place

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
2.2

Secure

Confirm SLEEQ EXT is firmly secured to Wing on both sides.

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
2.3

Re-attach

If included, re-Attach AeroFlex Anterior Panel. This support panel may be moved left or right to ensure it is centered on patient's abdomen.

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SLEEQ-Anwendung

Application de la ceinture SLEEQ

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Refer to Instructions for Use included with every SLEEQ brace to apply as usual.

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LIMITED WARRANTY

Thuasne, LLC warrants to the original purchaser that this product is free from defects in materials and workmanship. Normal wear and tear during the use of a product is not considered a defect. This Limited Warranty is effective from the date of delivery to the original purchaser. This warranty does not apply if this product:

- (1) has been altered in any way;
 - (2) was not purchased from Thuasne, LLC or one of its authorized distributors;
 - (3) is not used in accordance with Thuasne, LLC's Instructions for Use (IFU).
- Under this Limited Warranty Thuasne, LLC's sole obligation shall be to repair or replace the defective product or associated part(s) at no charge. This Limited Warranty does not cover damage as a result of mishandling, accidents, improper use or alteration of the product. This warranty is in place of all other warranties, either express or implied, including but not limited to warranties of merchantability or fitness for a particular purpose. Thuasne, LLC shall not be liable under this or any implied warranty for any personal injury, property damage, or special, incidental, or consequential damages arising out of, related, or incidental to the use of the product even if Thuasne, LLC has been informed of the prospect of such possible loss or damage, unless state law otherwise prohibits this exclusion. User should seek advice from their healthcare provider immediately if user encounters any kind of adverse reaction related to the use of this product. User should consult with their healthcare provider with any questions or concerns relating to the proper use of this product for their specific medical condition. IMPORTANT- Read all Instructions for Use (IFU) and any other associated materials to ensure proper use.

THUASNE LLC
4615 Shepard Street
Bakersfield, CA 93313
United States of America

Warnings

- Carefully read use/care instructions and warnings prior to use.
- If increased pain, swelling, sensation changes, or any adverse reactions are experienced while using this product, immediately consult your medical professional.
- For single patient use only.

Composition

Nylon, polyester, polyurethane

Care/Maintenance

Refer to the maintenance instructions of the SLEEQ brace.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

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MEDENVOY GLOBAL B.V.
THUASNE UK Ltd