


BOA SI

The BOA SI provides support for the Sacro-iliac joint. Excellent adjunct to post-injection pain following an S.I. procedure. Reduces pain and swelling for strained lax Sacro-iliac ligaments.

***L0621****BOA Low Profile 8"**


The BOA Low Profile brace provides rigid anterior and posterior compression. The 8" model will accommodate most patients.

***L0627****BOA Slim Panel XTR**

The BOA Back Brace with Slim Panel XTR limits extension and provides excellent support from T9-S1.

***L0631****BOA w/Chairback**

The BOA with Chairback combines the BOA Low Profile brace (pictured above), with a rigid chairback frame to provide spinal support from T9 - S1.

***L0637****BOA Duel TLSO**


The BOA Duel TLSO is a unique prefabricated orthosis for post-operative, post-traumatic, or chronic spinal pathology. The Duel TLSO provides full-circumferential rigid support from the scapular thoracic region (T2 / T23) of the spine down to the Sacrocygeal junction (S1). The BOA Duel restricts forward flexion and encourages a hyperextension posture. This posture increases lumbar lordosis and resists lateral and rotary movements of the spine.

Rigid Anterior & Posterior Panels

- Restrict forward flexion and trunk rotation
- Are easily removed and thermoformed for custom fit

Telescoping Sternal Y-Bar

- Can be bent to desired angle for optimal compression
- Easily adjusted to desired height, with lock screw for patient safety and compliance

***L0460**

General Information

Ordered by : _____

Phone: _____

Patient's Name: _____

BILLING: P.O. Number _____

Townsend Account Number: _____

Bill To: _____

Address: _____

City: _____

State: _____ Zip Code: _____ Country: _____

Phone: _____ Fax: _____

Ship To: _____

Address: _____

City: _____

State: _____ Zip Code: _____ Country: _____

Phone: _____ Attention: _____

Shipping Preference:

Ground 2-Day P.M. 2-Day A.M.

Next Day P.M. Next Day A.M.

(If no preference is indicated, this order will be shipped 2 Day P.M.)

Note: We do not ship new or repaired products directly to patients.

Sizing Reference

Model	Part	X-Small	Small	Medium	Large	X-Large	XX-Large	XXX-Large
BOA Slim Panel XTR LP (8")	51023-x	20" - 25"	25" - 30"	30" - 35"	35" - 40"	40" - 45"	45" - 50"	50" - 55"
BOA LP with Chairback (8")	51020-x	N/A	25" - 30"	30" - 35"	35" - 40"	40" - 45"	45" - 50"	50" - 55"
BOA LP (8")	51010-x	20" - 25"	25" - 30"	30" - 35"	35" - 40"	40" - 45"	45" - 50"	50" - 55"
BOA SI Belt	51024-x	N/A	25" - 35"	35" - 45"	45" - 55"	N/A	N/A	N/A
BOA Duel TLSO	51050-x	N/A	25" - 30"	30" - 35"	35" - 42"	42" - 48"	48" - 57"	N/A

BOA Spinal Braces

Models

BOA Slim Panel XTR LP (8") -- 51023--L0631
Quantity: XS ___ S ___ M ___ L ___ XL ___ XXL ___ XXXL ___

BOA LP with Chairback (8") -- 51020--L0637
Quantity: S ___ M ___ L ___ XL ___ XXL ___ XXXL ___

BOA LP (8") -- 51010--L0627
Quantity: XS ___ S ___ M ___ L ___ XL ___ XXL ___ XXXL ___

BOA SI Belt -- 51024--L0621
Quantity: S ___ M ___ L ___

BOA Duel TLSO -- 51050--L0460
Quantity: S ___ M ___ L ___ XL ___ XXL ___

Accessories

Standard Chairback -- 51015
Quantity: S ___ M ___ L ___ XL ___

Short Chairback-- 51016
Quantity: M ___ L ___ XL ___ XXL ___ XXXL ___

8" Slim Panel XTR (panel only) -- 510230

7° Lordotic Insert (pair) -- 510150

10° Lordotic Insert (pair) -- 510151

12° Lordotic Insert (pair) -- 510152



Townsend Design

4615 Shepard St., Bakersfield, CA, 93313
Phone: 800.432.3466 or 661.837.1795; Fax: 800.798.2722

Ordering Instructions

Please complete and fax this form to 800.798.2722 (24-hours a day). If you are calling in your order, this form indicates the options and information that will be required by our staff. For phone orders, please call 800.700.2722 between 6:00 a.m. and 4:00 p.m. (PST).