



LOMBA[®] SPINAL BRACES

PDAC Approved Lumbar Orthoses



THUASNE

LOMBASKIN PRO® LOMBATECH PRO®


PREMIUM WOVEN
COMPRESSION BRACES

Anatomically Adaptive Technology

- Woven fabric conforms perfectly to the anatomy to optimize fit, comfort, and migration resistance
- Complete modularity allows adapting from maximum stabilization to full step-down for progressive rehabilitation

Combitex® Patented One-way Stretch Weave

- Horizontal stretch with vertical inelasticity focuses compression for enhanced stability



Exclusive one-way
stretch weave
maximizes stability

Comfort and Compliance

- Breathable fabrics with seamless edges reduce moisture and skin pressure
- Visually appealing, ultra-low profile design



Removable Rigid
Back Panel

LombaSkin Pro®

PDAC approved L0627/L0642, L0631/L0648

- 3 sizes accommodate waists from 26-62 inches and minimize inventory
- Removable rigid back panels insert into posterior pocket
- Dual elastic straps multiply compressive forces



Lateral
Panels

Additional
Compression
Strap Included

LombaTech Pro®

PDAC approved L0631/L0648, L0637/L0650

- Secondary strap with user friendly hand loops provides added compression
- Removable anterior, posterior and lateral panels
- Adjustable posterior stays for enhanced support



4615 Shepard Street Bakersfield, CA 93313
Telephone 800.432.3466 - Fax 800.798.2722
www.ThuasneUSA.com