

Sport Series

Ligament Bracing Solutions

Routinely Prescribed For...

Patients who have experienced mild ligament injuries and/or reconstructive knee surgery.

This Product Can Be Ordered...

Offered With These Standard Features...

- Compression Sleeve
- Patented Townsend Motion Hinges
- Rigid Aluminum Bars
- Reinforced Nylon Hinge Pockets
- Condylar Pads

Available Options Include...

- Latex Free Neoprene Sleeve Material
- Fiber Cool Sleeve Material (Anterior Closure Models)
- Pull On Or Anterior Closure Models
- 12 Inch (2 Strap) and 16 Inch (4 Strap) Models
- Two Hinge Options
- Extension Stop Kit
- Flexion Stop Kit
- Patella Buttress (Pull-On Sleeves)



Townsend's Sport Series braces are designed for maximum durability and achieve more substantial control than traditional soft shell products. The patented low profile hinges track with the knee and are fully contained inside reinforced pockets so there is no exposed metal.

Hinge & Sleeve Options

Sport Series braces are available in a variety of configurations. The optional Fiber Cool sleeves feature an open cell material that reduces sweating and enhances comfort. The anterior closure models are easier to put on and allow for individual tightening of the thigh and tibia sections. The longer four strap models are recommended for taller patients and additional stabilization. The standard hinges have a fixed nonadjustable 5 degree extension stop; the TM5+ Hinges include adjustable extension stops, and optional flexion stops.



TM5+

Standard

 **Townsend**

www.townsenddesign.com

Toll Free: 800.432.3466 Ph: 661.837.1795 Fax: 661.837.0613

General Information

Ordered by : _____

Phone: _____

Patient's Name: _____

BILLING: P.O. Number _____

Townsend Account Number: _____

Bill To: _____

Address: _____

City: _____

State: _____ Zip Code: _____ Country: _____

Phone: _____ Fax: _____

Ship To: _____

Address: _____

City: _____

State: _____ Zip Code: _____ Country: _____

Phone: _____ Attention: _____

Fit Date: If known, please indicate the date you are scheduled to fit the patient: _____

Shipping Preference:

Ground 2-Day P.M. 2-Day A.M.

Next Day P.M. Next Day A.M.

(If no preference is indicated, this order will be shipped 2 Day P.M.)

Note: We do not ship new or repaired products directly to patients.

Sizing Reference

With the patient's leg at full extension, measure circumference around the center of the affected knee.

Small: 13" to 14" circumference

Medium: 14" to 15" circumference

Large: 15" to 17" circumference

XL: 17" to 19" circumference

XXL: 19" to 21" circumference

XXXL: 21" to 22" circumference

XXXXL: 22" to 23" circumference (neoprene sleeves only)

Ordering Instructions

Please complete and fax this form to 800.798.2722 (24-hours a day). If you are calling in your order, this form indicates the options and information that will be required by our staff. For phone orders, please call 800.700.2722 between 6:00 a.m. and 4:00 p.m. (PST).

Townsend Sport Series

*Indicates additional charges apply

Hinge Descriptions

Standard nonadjustable hinges with a 5 degree built-in extension stop

TM5+ Hinges*

Includes Extension Stop Kit: 0, 5, 10, 15 & 20 degrees

Optional Flexion Stop Kit: 0, 30, 45, 60, 75, 90 and 110 degrees



Standard



TM5+

Fiber Cool Sleeve (breathable non-neoprene material)

Sport XII-AC (12" anterior closure sleeve, 2 straps)

Size/Quantity: S ___ M ___ L ___ XL ___ XXL ___ XXXL ___

Hinge Choice: 5 Degree Hinges TM5+ Hinges*

Sport XVI-AC (16" anterior closure sleeve, 4 straps)

Size/Quantity: S ___ M ___ L ___ XL ___ XXL ___ XXXL ___

Hinge Choice: 5 Degree Hinges TM5+ Hinges*

Neoprene Sleeve

Sport XII-AC (12" anterior closure sleeve, 2 straps)

Size/Quantity: S ___ M ___ L ___ XL ___ XXL ___ XXXL ___ XXXXL ___

Hinge Choice: 5 Degree Hinges TM5+ Hinges*

Sport XII-PO (12" pull on sleeve, 2 straps)

Size/Quantity: S ___ M ___ L ___ XL ___ XXL ___ XXXL ___ XXXXL ___

Hinge Choice: 5 Degree Hinges TM5+ Hinges*

Sport XVI-AC (16" anterior closure sleeve, 4 straps)

Size/Quantity: S ___ M ___ L ___ XL ___ XXL ___ XXXL ___ XXXXL ___

Hinge Choice: 5 Degree Hinges TM5+ Hinges*

Sport XVI-PO (16" pull on sleeve, 4 straps)

Size/Quantity: S ___ M ___ L ___ XL ___ XXL ___ XXXL ___ XXXXL ___

Hinge Choice: 5 Degree Hinges TM5+ Hinges*

Parts/Accessories

Replacement Sleeve (indicate size and model):

Townsend Design

4615 Shepard St., Bakersfield, CA, 93313

Phone: 800.432.3466 or 661.837.1795; Fax: 800.798.2722

Special Instructions: _____
