

LOMBA SPINAL BRACES

Premium Woven Compression Orthoses



LombaSkin Pro®

L0627/L0642, L0631/L0648

- Removable rigid back panels insert into posterior pocket
- 3 sizes accommodate waists from 26" - 62" to minimize inventory
- Dual elastic straps multiply compressive forces

LombaTech Pro®

L0627/L0642, L0631/L0648, L0637/L0650

- Removable anterior, posterior and lateral panels
- Secondary strap with user friendly hand loops provides added compression
- Adjustable posterior stays for enhanced support

Lomba Model	Part #	X-Small (X = 1)	Small (X = 2)	Medium (X = 3)	Large (X = 4)	X-Large (X = 5)
LombaSkin Pro (AP)	51030-x	N/A	26" - 36"	36" - 47"	48" - 62"	N/A
LombaSkin Pro (AP+)	51032-x	N/A	26" - 36"	36" - 47"	48" - 62"	N/A
LombaTech Pro (AP)	51060-x	23½" - 30"	30½" - 37"	37½" - 45"	45½" - 55"	55½" - 63"
LombaTech Pro (AP+)	51061-x	23½" - 30"	30½" - 37"	37½" - 45"	45½" - 55"	55½" - 63"
LombaTech Pro (APL)	51062-x	23½" - 30"	30½" - 37"	37½" - 45"	45½" - 55"	55½" - 63"

*Measure Waist Circumference

Dorso Rigid

Premium Hyperextension Orthosis

L0472

- Optimized fit and control
- Lightweight telescoping aluminum alloy frame
- Back pad positions horizontally or vertically
- Articulating sternal /pelvic pads

Part # 089001



	Pelvic Circumference	Pubic Symphysis to Sternum Height
S	23½" - 29½"	12½"
MS	29½" - 35½"	12½"
M	29½" - 35½"	13¾"
LS	35½" - 41¼"	13¾"
L	35½" - 41¼"	15"
XL	41¼" - 45¼"	16¼"

LombaMum

Maternity Lumbar Belt

- For low back, pelvic and sacro-iliac pain
- Strapping system adjusts amount of support
- High quality materials for breathability and comfort

Part # 080501

Universal sizing: (33 1/2" - 49 1/4")



CERVICAL BRACING



The Evolution of the Cervical Collar

Eclipse® Universal Cervical Collar L0172, L0174

- Effortless chin height adjustment for easy fitting and maximum compliance
- Fifteen height settings for precise positioning
- Quick-Snap fitting disks ensure proper patient application
- Anatomically molded, structurally reinforced frame for superior strength and fit

Eclipse (Standard) Part # 910010

Eclipse EXT (Extended Thoracic) Part # 910020

Universal sizing
(Wide back panel for neck circumference greater than 19")
Replacement pad sets available

Ortel C42 Vario® L0174

- Adjustable dual shell design provides immobilization of the cervical spine
- Dual height adjustment knobs range from 3½" to 6 ¼"
- X-ray friendly, radio-transparent design with large tracheal opening
- High quality, breathable, removable liner

Part # T492801 Universal Sizing



Thuasne C1

- Soft foam support for the cervical spine
- High-quality foam, removable and washable cover
- Recommended for short-term use

Part # 135002

Size	Collar Height
1	2¾"
2	3½"
3	4¾"



Ortel C1 Junior®

- Pediatric soft foam support for the cervical spine
- Removable rigid insert for adjustability

Part # 2620001004

Universal Sizing:
(9½" - 10½") x 2½" high



Ortel C2 Plus

- Semi-rigid internal structure
- Anatomical shape with rounded edges
- Removable and washable cover

Part # 239203

Available in three heights:
3", 3½", and 4¼"

Size	Neck Circumference
1	11" - 13"
2	13½" - 15¼"
3	15½" - 18"



SPINAL BRACING PRODUCTS



4615 Shepard St. • Bakersfield, CA 93313
800.432.3466 toll free
800.798.2722 toll free fax




1000 Calle Cordillera • San Clemente, CA 92673
855.784.6600 toll free
855.784.6601 toll free fax

www.ThuasneUSA.com

© 2016 THUASNE USA All rights reserved.



SLEEQ Spinal Therapy System



The Most Comfortable & Easy-to-Adjust Spinal Compression Braces

Universal Sizing

Extendable panels accommodate waist sizes from 26" to 62" – wings can be angled to optimize fit and comfort for all body types.

Dual-Zone Technology

Delivers optimal non-elastic compression with relief at the ribs and hips for unparalleled comfort and compliance – even while seated.

Partner Solutions

SLEEQ is supported by industry leading marketing, reimbursement, and clinical education programs.

Proudly Made in USA



SLEEQ's ingenious universally adjustable design reduces inventory and provides exceptional compression, pain relief and comfort for maximum patient compliance.



SLEEQ AP
L0627/L0642
Low-profile sagittal (anterior-posterior) control lumbar orthosis
Part # 901100



SLEEQ AP+
L0631/L0648
Enhanced profile sagittal (anterior-posterior) control lumbo-sacral orthosis
Part # 901200



SLEEQ APL
L0637/L0650
Enhanced profile sagittal-coronal (anterior-posterior-lateral) control lumbo-sacral orthosis
Part # 901300



SLEEQ FLEX
L0456/L0457
Flexible sagittal (anterior-posterior) control thoraco-lumbo-sacral orthosis
Part # 901400

SLEEQ Accessories

SLEEQ EXT

Extension panel lengthens any SLEEQ to fit waists from 54"– 62"
Part # 901070



Q^othermo

Premium cold/heat pack designed to be used with all SLEEQ braces for multi-modality pain therapy
Part # 950200



(((SLEEQ))) compatible



BOA braces provide excellent semi-rigid support via mechanical advantage pulley systems and medical grade plastics in well proven traditional brace designs.



BOA SI
L0621
Ultra low profile support for the sacro-iliac joint available in three sizes
Part # 51024-X



BOA LP 8"
L0627/L0642
Low profile anterior-posterior control
Part # 51010-X



BOA Slim Panel XTR
L0631/L0648
Semi-rigid anterior-posterior control
Part # 51023-X

BOA w/Chairback

Semi-rigid anterior-posterior with lateral control
Part # 51020-X



BOA Duel TLSO

Full-circumferential rigid support for the thoraco-lumbar spine
Part # 51050-X



BOA Model	Part #	X-Small (X = 2)	Small (X = 3)	Medium (X = 4)	Large (X = 5)	X-Large (X = 6)	XX-Large (X = 7)	XXX-Large (X = 8)
BOA Duel TLSO	51050-x	N/A	25"- 30"	30"- 35"	35"- 42"	42"- 48"	48"- 57"	N/A
BOA LP with Chairback (8")	51020-x	20"- 25"	25"- 30"	30"- 35"	35"- 40"	40"- 45"	45"- 50"	50"- 55"
BOA LP (8")	51010-x	20"- 25"	25"- 30"	30"- 35"	35"- 40"	40"- 45"	45"- 50"	50"- 55"
BOA SI Belt	51024-x	N/A	25"- 35"	35"- 45"	45"- 55"	N/A	N/A	N/A
BOA Slim Panel XTR (8")	51023-x	20"- 25"	25"- 30"	30"- 35"	35"- 40"	40"- 45"	45"- 50"	50"- 55"

*Measure Waist Circumference

BOA Accessories	Part #	Small (X = 4)	Medium (X = 5)	Large (X = 6)	X-Large (X = 7)
Standard Chairback	51015-x	25"- 30" /13"	30"- 35" /14"	35"- 40" /14½"	40"- 55" /14½"
Short Chairback	51016-x	N/A	30"- 35" /13"	35"- 40" /13"	40"- 55" /13"

The Short Chairback is designed for shorter torsos with a 13" frame height for all sizes.

8" Slim Panel XTR (panel only)	510230	Panel can be added to BOA LP	7°, 10°, 12° Lordotic inserts available.
--------------------------------	--------	------------------------------	--