



**THUASNE**



## **DORSO RIGID<sup>®</sup> 35**

PDAC Approved L0472

Control thoracolumbar flexion with the Dorso Rigid 35 premium hyperextension brace.



**TOWNSEND**  
THUASNE USA



**QUINN MEDICAL**  
THUASNE USA

# Dorso Rigid 35

Premium Hyperextension Orthosis



## Easily Optimize Fit And Control

- Adjustable telescoping aluminum alloy frame
- Thoracic band widens without removing screws, and adjusts angle in 7.5° increments
- Malleable frame for optimized fit
- Back pad positions horizontally or vertically



## ARTICULATING PELVIC BAND AND STERNAL PAD

- Spring assisted pelvic band articulates or locks at 0° or 18°
- Sternal pad self-adjusts for comfortable control

## EASY TO USE LEVER LOCK

- Effortless application and tightening for enhanced compliance
- Safe and secure latch closure
- Quick connect strap button with hand loop for easy application



Size	Pelvis Circumference	Pubic Symphysis To Sternum Height	
SM	23 ½ – 29 ½ inch	12 ½ inch	
MS	29 ½ – 35 ½ inch	12 ½ inch	
MD	29 ½ – 35 ½ inch	13 ¾ inch	
LS	35 ½ – 41 ¼ inch	13 ¾ inch	
LG	35 ½ – 41 ¼ inch	15 inch	
XL	41 ¼ – 45 ¼ inch	16 ¼ inch	



## Water Resistant Materials

- Brace can be immersed in water for showering or pool rehabilitation
- Stainless steel screws, aluminum frame, and plastic hardware
- Water repellent pads