

Ordered by: _____ Phone # (_____) _____

BILLING: P.O. Number _____ Townsend Account # _____

Bill To: _____ Ship To: _____

Address: _____ Address: _____

City: _____ City: _____

State: _____ Zip Code: _____ Country: _____ State: _____ Zip Code: _____ Country: _____

Phone: (_____) _____ Fax: (_____) _____ Phone: (_____) _____ Fax: (_____) _____

Shipping Preference: Ground 2-Day P.M. 2-Day A.M. Next Day P.M. Next Day A.M.

(If no preference is indicated, this order will be shipped 2 Day P.M.) Note: We do not ship products directly to patients.

Patient's Last Name: _____

Patient's First Name: _____

Measuring Information

Measure arm from elbow crease to base of the index finger:
S/M up to 14 inches; L/XL up to 15 ½ inches.

**Low Profile ROM Elbow (Dual Hinges) With Shells And
Optional Hand Extension**

- Flexion can be limited at 90, 110 and 120 degrees, and extension can be stopped at 0, 15, 30, 45 and 90 degrees.
- The screw stops are adjusted using a small flathead screwdriver.
- The hinge can be locked at 90 degrees.

Quantity: ___ S/M ___ L/XL

Install Optional Hand Extension? ___ Yes ___ No

Include Shoulder Sling Strap?* ___ Yes ___ No

Multi-Adjustment ROM Elbow

- Comes with four straps and individual black pads that fit under each strap.
- Tool-less flexion and extension stops are in 10 degree increments.
- Locking is possible by offsetting the flexion and extension stops at the same degree.

Single Hinge (SH) Model

Quantity: ___ LEFT ARM Length: ___ 15 inch ___ 16 inch

Quantity: ___ RIGHT ARM Length: ___ 15 inch ___ 16 inch

Dual Hinge (DH) Model

Length (Indicate quantity): ___ 15 inch ___ 16 inch

ROM-LBO

- The LBO comes with dual or single uprights, four straps, tool less adjustability, and 8 individual Breeze pads.
- Tool-less flexion and extension stops are in 10 degree increments.
- The LBO can be locked at 90 degrees for post-surgical protocols.
- Adjustable length 14 ½" – 17"

Single Hinge (SH) Model

Quantity: ___ LEFT ARM

Quantity: ___ RIGHT ARM

Dual Hinge (DH) Model

Quantity: ___ Universal

Please complete and fax this form to 800.798.2722 (24-hours a day). If you are calling in your order, this form indicates the options and information that will be required by our staff. For phone orders, please call 800.700.2722 between 6:00 a.m. and 4:00 p.m. (PST).

*Indicates additional charges apply