

Ordered by: \_\_\_\_\_ Phone # ( \_\_\_\_\_ ) \_\_\_\_\_

BILLING: P.O. Number \_\_\_\_\_ Townsend Account # \_\_\_\_\_

Bill To: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_

State: \_\_\_\_\_ Zip Code: \_\_\_\_\_ Country: \_\_\_\_\_

Phone: ( \_\_\_\_\_ ) \_\_\_\_\_ Fax: ( \_\_\_\_\_ ) \_\_\_\_\_

Ship To: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_

State: \_\_\_\_\_ Zip Code: \_\_\_\_\_ Country: \_\_\_\_\_

Phone: ( \_\_\_\_\_ ) \_\_\_\_\_ Fax: ( \_\_\_\_\_ ) \_\_\_\_\_

Shipping Preference:  Ground  2-Day P.M.  2-Day A.M.  Next Day P.M.  Next Day A.M.  
 (If no preference is indicated, this order will be shipped 2 Day P.M.) Note: We do not ship products directly to patients.

Patient's Last Name: \_\_\_\_\_

Patient's First Name: \_\_\_\_\_

Male  Female Age \_\_\_\_\_

Weight \_\_\_\_\_ (LBS) Height \_\_\_\_\_ (IN) Leg:  Left  Right

**Compartment**

- Medial Compartment  Universal  
(select when ordering stock inventory)  
 Medial OA pads installed, Lateral OA pads included in kit
- Lateral Compartment

**Shell Length** Thigh shell extends 7" to 9" on all braces

**Tibia**  6 Inch  7 Inch  8 Inch

**TM5+ Hinges** Includes Extension Stops

**Optional Flexion Stop Kit\***

**Standard Colors** in stock

- Argento Gray  Textured Black

**Non-stock Colors** Brace may ship next business day

- Bengal Silver  Electric Blue  Sparkle Copper  
 Bengal Yellow  Sky Blue  Sparkle Red  
 Bengal White  Violet  
 Gloss Black  Candy Green

**PCL Strap\*** (combined instabilities) Or  **CS Package\***

**Anti-Migration Silicon Infused Strap Pads\***

**Spooner Patella Stabilizing Attachment\***

**Easy Release Buckles\***  
 (Not available on the distal thigh strap when ordered with CS Package)

**Brace Cover\***

- Posterior Closure  Pull On

**Undersleeves\***

- 18" Cotton  18" Neoprene  22" Neoprene

**Thigh Sleeves\***

1/16 Comfort Thigh Sleeve

**Size**


- X-Small  Small  Medium  
 Large  X-Large  XX-Large

**Sizing Reference**

	6" above	M-L width	6" below
X-Small	NOT AVAILABLE	3" to 3.5"	11" to 12.5"
Small	15.5" to 18.5"	3.5" to 4"	12.25" to 13.75"
Medium	18.5" to 21"	4" to 4.5"	13.25" to 15"
Large	21" to 23.5"	4.5" to 5"	14.25" to 15.75"
X-Large	23.5" to 25"	5" to 5.5"	15" to 17"
XX-Large	25" to 28"	5.5" to 6"	17" to 19"

**3 Measurements: "Customized" Assembly** (No Added Charge)

If your patient has proportional leg sizing (see sizing reference, above) or if you are ordering a brace for stock inventory please select from the size options. However, if you would prefer to have Townsend customize the assembly of your patient's brace at no addition charge please provide leg measurements beside the illustration, below.



Please complete and fax this form to 800.798.2722 (24-hours a day). If you are calling in your order, this form indicates the options and information that will be required by our staff. For phone orders, please call 800.700.2722 between 6:00 a.m. and 4:00 p.m. (PST).

\*Indicates additional charges apply