




The Air Townsend is a traditional full shell design that provides maximum control and protection for ligament instabilities.



The Air Lite features an adjustable anterior tibia strap and narrower bands to reduce the weight and profile of the brace.



## Optional Shell Configuration

Optional single upright KAFO with leg extension and heel cup for rotary control or to eliminate any potential for brace migration.

# Air Townsend & Air Lite

## Ligament Bracing Solutions

### Routinely Prescribed For...

Patients who desire or require maximum control after a ligament injury and/or reconstructive knee surgery.

### This Product Can Be Ordered...

- Custom from a cast, digital scan or KneePix

### Offered With These Standard Features...

- Custom Molded Laminated Carbon Graphite Shells
- Patented Townsend Motion TM5+ Hinges
- Extension Stop Kit
- Patented Synergistic Suspension Strap
- Anti-Migration Strap Padding
- Anti-Rotation Tibia Shell Bolster (Air Townsend)

### Available Options Include...

- ACL Or Combined Instabilities
- Virtually Any Shell Length Over 13 Inches
- Hinge Material -- Aluminum Or Stainless
- Extension Assist
- Special Trim Lines For Skiers, Rodeo, Hockey Goalies, And Full Figure (Air Townsend)

[www.townsenddesign.com](http://www.townsenddesign.com)

Toll Free: 800.432.3466 Ph: 661.837.1795 Fax: 661.837.0613

 **Townsend**