

Ordered by: \_\_\_\_\_ Phone # ( \_\_\_\_\_ ) \_\_\_\_\_

BILLING: P.O. Number \_\_\_\_\_ Townsend Account # \_\_\_\_\_

Bill To: \_\_\_\_\_

Ship To: \_\_\_\_\_

Address: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_

City: \_\_\_\_\_

State: \_\_\_\_\_ Zip Code: \_\_\_\_\_ Country: \_\_\_\_\_

State: \_\_\_\_\_ Zip Code: \_\_\_\_\_ Country: \_\_\_\_\_

Phone: ( \_\_\_\_\_ ) \_\_\_\_\_ Fax: ( \_\_\_\_\_ ) \_\_\_\_\_

Phone: ( \_\_\_\_\_ ) \_\_\_\_\_ Fax: ( \_\_\_\_\_ ) \_\_\_\_\_

Shipping Preference:  Ground  2-Day P.M.  2-Day A.M.  Next Day P.M.  Next Day A.M.  
 (If no preference is indicated, this order will be shipped 2 Day P.M.) Note: We do not ship products directly to patients.

**Patient's Last Name:** \_\_\_\_\_

**Patient's First Name:** \_\_\_\_\_

Male  Female Age \_\_\_\_\_ Height \_\_\_\_\_ Weight \_\_\_\_\_

**Leg**  Left  Right

**TM5+ Hinges** (Includes Extension Stops)

Optional Flexion Stop Kit\*

**Options**

Spooner Patella Attachment\*

Motocross Knee Guard\*

**Anti-Migration Silicon Infused Strap Pads\***

**Brace Cover\***

Posterior Closure

Pull On

**Undersleeves\***

18" Cotton  18" Neoprene  22" Neoprene

**Thigh Sleeves\***

1/16 Comfort Thigh Sleeve


**Size**

Small  Medium  Large  X-Large

**Sizing Reference**

	6" above	M-L width	6" below
<b>Small</b>	15.5" to 18.5"	3.5" to 4"	12.25" to 13.75"
<b>Medium</b>	18.5" to 21"	4" to 4.5"	13.25" to 15"
<b>Large</b>	21" to 23.5"	4.5" to 5"	14.25" to 15.75"
<b>X-Large</b>	23.5" to 25"	5" to 5.5"	15" to 17"

Please provide the three measurements, below, so we can determine the most appropriate standard size BOLD brace to ship. The knee width M-L is the PRIMARY measurement because the central aluminum structure is very rigid and not easily adjusted. The leg circumferences are SECONDARY reference measurements, as the semi-rigid nylon shells will generally contour to the leg when the thigh and tibia straps are tightened. For a disproportional leg, consider using a "custom assembled" Townsend Rebel or Rebel Pro brace.



Please complete and fax this form to 800.798.2722 (24-hours a day). If you are calling in your order, this form indicates the options and information that will be required by our staff. For phone orders, please call 800.700.2722 between 6:00 a.m. and 4:00 p.m. (PST).

\* Indicates additional charges apply