

Account Contact Information

Name: _____ Email: _____ Phone: _____

Billing and Shipping

PO# _____ Billing Account #: _____ Shipping Account #: _____

Billing Address: _____ Shipping Address: _____

- Shipping Preference**
- Ground
 - Next Day A.M.
 - Next Day P.M.
 - 2-Day A.M.
 - 2-Day P.M.

City: _____ City: _____

State _____ Zip _____ State: _____ Zip: _____

(If no preference is indicated, this order will be shipped 2 Day P.M.) Note: We do not ship products directly to patients.

Patient Information

Patient's Last Name: _____

Patient's First Name: _____

Measuring Information

Measure arm from elbow crease to base of the index finger:
S/M up to 14 inches; L/XL up to 15 ½ inches.

Low Profile ROM Elbow (Dual Hinges) With Shells And Optional Hand Extension

- Flexion can be limited at 90, 110 and 120 degrees, and extension can be stopped at 0, 15, 30, 45 and 90 degrees.
- The screw stops are adjusted using a small flathead screwdriver.
- The hinge can be locked at 90 degrees.

Quantity: _____ S/M _____ L/XL

Install Optional Hand Extension? _____ Yes _____ No

Include Shoulder Sling Strap? * _____ Yes _____ No

Multi-Adjustment ROM Elbow

- Comes with four straps and individual black pads that fit under each strap.
- Tool-less flexion and extension stops are in 10 degree increments.
- Locking is possible by offsetting the flexion and extension stops at the same degree.

Single Hinge (SH) Model

Quantity: _____ LEFT ARM Length: _____ 15 inch _____ 16 inch

Quantity: _____ RIGHT ARM Length: _____ 15 inch _____ 16 inch

Dual Hinge (DH) Model

Length (Indicate quantity): _____ 15 inch _____ 16 inch

ROM-LBO

- The LBO comes with dual or single uprights, four straps, tool less adjustability, and 8 individual Breeze pads.
- Tool-less flexion and extension stops are in 10 degree increments.
- The LBO can be locked at 90 degrees for post-surgical protocols.
- Adjustable length 14 ½" – 17"

Single Hinge (SH) Model

Quantity: _____ LEFT ARM

Quantity: _____ RIGHT ARM

Dual Hinge (DH) Model

Quantity: _____ Universal

*Indicates additional charges apply