

Account Contact Information

Name: _____ Email: _____ Phone: _____

Billing and Shipping

PO# _____ Billing Account #: _____ Shipping Account #: _____

Billing Address: _____ Shipping Address: _____

City: _____ City: _____

State: _____ Zip: _____ State: _____ Zip: _____

Shipping Preference

- Ground
- Next Day A.M.
- Next Day P.M.
- 2-Day A.M.
- 2-Day P.M.

(If no preference is indicated, this order will be shipped 2 Day P.M.) Note: We do not ship products directly to patients.

Patient Information

Patient's Last Name: _____

Patient's First Name: _____

Fiber Cool Sleeve *(breathable non-neoprene material)*

Sport XII-AC (12" anterior closure sleeve, 2 straps)

Hinge Choice: 5 Degree Hinges TM5+ Hinges*
 Size/Quantity: ___ S ___ M ___ L ___ XL
 ___ 2XL ___ 3XL ___ 3XL

Sport XVI-AC (16" anterior closure sleeve, 4 straps)

Hinge Choice: 5 Degree Hinges TM5+ Hinges*
 Size/Quantity: ___ S ___ M ___ L ___ XL
 ___ 2XL ___ 3XL ___ 3XL

Neoprene Sleeve

Sport XII-AC (12" anterior closure sleeve, 2 straps)

Hinge Choice: 5 Degree Hinges TM5+ Hinges*
 Size/Quantity: ___ S ___ M ___ L ___ XL
 ___ 2XL ___ 3XL ___ 3XL ___ 4XL

Sport XVI-AC (16" anterior closure sleeve, 4 straps)

Hinge Choice: 5 Degree Hinges TM5+ Hinges*
 Size/Quantity: ___ S ___ M ___ L ___ XL
 ___ 2XL ___ 3XL ___ 3XL ___ 4XL

Hinge Descriptions

Standard Nonadjustable Hinges

5 degree built-in extension stop

TM5+ Hinges*

Includes Extension Stop Kit: 0, 5, 10, 15 & 20 degrees

Sizing Reference		<i>With the patient's leg at full extension, measure circumference around the center of the affected knee.</i>
	circumference	
Small	13" to 14"	
Medium	14" to 15"	
Large	15" to 17"	
X-Large	17" to 19"	
2X-Large	19" to 21"	
3X-Large	21" to 22"	
4X-Large	22" to 23"	<i>(neoprene sleeves only)</i>

Replacement Sleeve *(indicate size and model):*

*Indicates additional charges apply