

**Account Contact Information**

Name: \_\_\_\_\_ Email: \_\_\_\_\_ Phone: \_\_\_\_\_

**Billing and Shipping**

PO# \_\_\_\_\_ Billing Account #: \_\_\_\_\_ Shipping Account #: \_\_\_\_\_

Billing Address: \_\_\_\_\_ Shipping Address: \_\_\_\_\_

**Shipping Preference**

- Ground
- Next Day A.M.
- Next Day P.M.
- 2-Day A.M.
- 2-Day P.M.

City: \_\_\_\_\_ City: \_\_\_\_\_  
 State: \_\_\_\_\_ Zip: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

*(If no preference is indicated, this order will be shipped 2 Day P.M.) Note: We do not ship products directly to patients.*

**Patient Information (If applicable)**

**Fit Date:** \_\_\_\_\_ **Patient Shoe Size:** \_\_\_\_\_

**Patient's Last Name:** \_\_\_\_\_

**Patient's First Name:** \_\_\_\_\_

**SpryStep® Original (weight limit = 250lbs)**

Right \_\_\_ XS \_\_\_ S \_\_\_ M \_\_\_ L \_\_\_ XL  
 Left \_\_\_ XS \_\_\_ S \_\_\_ M \_\_\_ L \_\_\_ XL

**SpryStep® Flex (weight limit = 250lbs)**

Right \_\_\_ XS \_\_\_ S \_\_\_ M \_\_\_ L \_\_\_ XL  
 Left \_\_\_ XS \_\_\_ S \_\_\_ M \_\_\_ L \_\_\_ XL

**SpryStep® Max (weight limit = 250lbs)**

Right \_\_\_ XS \_\_\_ S \_\_\_ M \_\_\_ L \_\_\_ XL  
 Left \_\_\_ XS \_\_\_ S \_\_\_ M \_\_\_ L \_\_\_ XL

**SpryStep® Plus (weight limit = 250lbs)**

Right \_\_\_ XS \_\_\_ S \_\_\_ M \_\_\_ L \_\_\_ XL  
 Left \_\_\_ XS \_\_\_ S \_\_\_ M \_\_\_ L \_\_\_ XL

**SpryStep® Pediatric (weight limit = 130lbs)**


Right \_\_\_ XS \_\_\_ S \_\_\_ M \_\_\_ L \_\_\_ XL  
 Left \_\_\_ XS \_\_\_ S \_\_\_ M \_\_\_ L \_\_\_ XL

**Accessories**

Refurbishment kit (pad, strap, alligator tab)\*

1. The forefoot section can be trimmed with scissors: Approximately 1 inch can be taken off the toe end, and ½ to ¼ inch can be trimmed from EACH side of the forefoot. The material that can be cut off the forefoot is tinted blue.
2. This indicates the shorter the footplate can be trimmed with scissors
3. This indicates the lowest the calf cuff can be trimmed with scissors

**A brace that has been trimmed cannot be returned for credit or replacement**

SpryStep® & SpryStep® Flex			Original	Flex
Size	Height	Foot Plate Length <sup>1</sup>		
XS	11 ¾"	8 ½ — 9 ½"		
SM	12 ⅝"	9 ⅞ — 10 ¼"		
MD	13 ⅜"	9 ⅝ — 10 ⅞"		
LG	14 ⅙"	10 ¼ — 11 ⅜"		
XL	15"	11 ⅜ — 12"		

SpryStep® Max		
Size	Height	Foot Plate Length <sup>1</sup>
XS	13 ¼"	7 ¾ — 8 ½"
SM	14 ½"	8 ¼ — 9"
MD	15 ½"	9 — 10"
LG	16 ½"	10 — 10 ¾"
XL	16 ½"	10 ¾ — 11 ½"

SpryStep® Plus		Calf circumference 1" below fibula head	
Size	Height to Tib Tuberosity	Foot Plate Length <sup>1</sup>	
XS	13"	8 ½ — 9 ⅝"	11 ½ — 14 ½"
S	14 ¼"	9 ⅞ — 10 ¼"	12 ¾ — 15 ¾"
M	15 ½"	9 ⅝ — 10 ⅞"	13 ½ — 17 ½"
L	16 ¾"	10 ¼ — 11 ⅜"	14 ¾ — 18 ¾"
XL	18"	11 ⅜ — 12"	15 ½ — 20 ½"

SpryStep® Pediatric				
Size	Foot Length <sup>2</sup>		Calf Height <sup>3</sup>	
	Min	Max	Min	Max
XS	3 ¾"	5 ¾"	6 ¼"	7"
SM	4 ¼"	6 ¼"	7"	7 ¾"
MD	5"	7"	7 ¾"	8 ½"
LG	5 ½"	7 ½"	8 ¾"	9 ½"
XL	6"	8"	9 ½"	10 ¼"

*Calf cuff depth = 2 ½" (the width of the calf cuff can be trimmed up to 1 ½" on each side)*

\*Indicates additional charges apply