

**Account Contact Information**

Name: \_\_\_\_\_ Email: \_\_\_\_\_ Phone: \_\_\_\_\_

**Billing and Shipping**

PO# \_\_\_\_\_ Billing Account #: \_\_\_\_\_ Shipping Account #: \_\_\_\_\_

**Shipping Preference**

Ground

Next Day A.M.

Next Day P.M.

2-Day A.M.

2-Day P.M.

Billing Address: \_\_\_\_\_ Shipping Address: \_\_\_\_\_

City: \_\_\_\_\_ City: \_\_\_\_\_

State \_\_\_\_\_ Zip \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

*(If no preference is indicated, this order will be shipped 2 Day P.M.) Note: We do not ship products directly to patients.*

**Patient Information**

**Fit Date:** \_\_\_\_\_

**Patient's Last Name:** \_\_\_\_\_

**Patient's First Name:** \_\_\_\_\_

Male  Female Age \_\_\_\_\_

Weight \_\_\_\_\_ (LBS) Height \_\_\_\_\_ (IN)

Leg:  Left  Right

Ligament:  ACL  PCL  LCL  MCL

Meniscus Damage:  Medial  Lateral

Surgeries (type/date): \_\_\_\_\_

AirTownsend  Air Lite

Model:  ACL  CombinedInstabilities  
 Thigh Strap\*  GraphiteBand\*

**Thigh Shell Length**

7"  8"  Other\_\_\_\_\_"

**Thigh Band Width**

1.5"  2"  2.5"  Other\_\_\_\_\_"

**Tibia Shell Length**

7"  8"  9"  Other\_\_\_\_\_"

Single Strut KAFO With Heel Cup  
*(Must Complete Additional Form For Brace Extension)*

**Hinges**

**TM5+ Hinges** — Includes extension stop kit  
 Optional flexion stop kit\*  
 Add optional extension assist bands/posts\*

**Hinge Material**

6061 Aluminum *(standard, if no hinge material is indicated)*  
 Stainless Steel\*

**Original Hinges\*** *Stainless (standard) includes extension stop kit*  
 Condylar Pads:  No  Bi-Lateral  Medial  Lateral  
 Optional flexion stop kit\*  
 Add optional extension assist bands/posts\*

**Tibia Shell Anti-Rotation Bolster (Air Townsend only)**

Standard  Minimal  Soft  No Bolster

**Special Trim Lines (Air Townsend only)**

Full Figure: For additional soft tissue containment  
 No tibia shell hole *(Townsend Original)*  
 Ski boot cut: 1/2 inch notched in distal tibia shell  
 Rodeo: Medial thigh cut-out for rider comfort  
 Customized Shell Design *(include instructions/drawing)*

**Synergistic Suspension Strap Attachment**

Standard attachment *(Lateral end recessed inside tibia shell)*  
 Lateral end attached to outer shell *(for prominent fibular head)*

**Color/Fabric Inlay**

Black  Royal Blue  Sheer Teal\*  
 Beige  Green  Sheer Purple\*  
 Gray  Burgundy  US Flag Fabric\*  
 Red  Clear Graphite  Fabric  
 Navy Blue  Sheer Red\* *-1 yard from patient\**

**Anti-Migration Silicon Enfused Strap Pads\***

**Spooner Patella Stabilizing Attachment\***

**Brace Cover\***

Pull On

**Undersleeves\***

18" Cotton  18" Neoprene  22" Neoprene

**Thigh Sleeves\***

1/8 Atrophy Thigh Sleeve  1/16 Comfort Thigh Sleeve

**Special Instructions:** \_\_\_\_\_

\*Indicates additional charges apply