

**Account Contact Information**

Name: \_\_\_\_\_ Email: \_\_\_\_\_ Phone: \_\_\_\_\_

**Billing and Shipping**

PO# \_\_\_\_\_ Billing Account #: \_\_\_\_\_ Shipping Account #: \_\_\_\_\_

**Shipping Preference**

Ground

Next Day A.M.

Next Day P.M.

2-Day A.M.

2-Day P.M.

Billing Address: \_\_\_\_\_ Shipping Address: \_\_\_\_\_

City: \_\_\_\_\_ City: \_\_\_\_\_

State \_\_\_\_\_ Zip \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

*(If no preference is indicated, this order will be shipped 2 Day P.M.) Note: We do not ship products directly to patients.*

**Patient Information**

**Fit Date:** \_\_\_\_\_

**Patient's Last Name:** \_\_\_\_\_

**Patient's First Name:** \_\_\_\_\_

Male  Female Age \_\_\_\_\_

**Weight** \_\_\_\_\_ (LBS) **Height** \_\_\_\_\_ (IN)

**Leg:**  Left  Right

**Patient's Diagnosis:** \_\_\_\_\_

**Surgeries (type/date):** \_\_\_\_\_

**Comments:** \_\_\_\_\_

**Model** *\*Indicates additional charges apply*

- Standard (with anterior thigh and calf shells)
- Polio (anterior thigh shell, posterior calf shell)

**Thigh Shell Length**

7"  8"  Other \_\_\_\_\_"

**Tibia Shell Length**

7"  8"  9"  Other \_\_\_\_\_"

Graphite Extension To Floor  
*(Must Complete Additional Form For Brace Extension)*

**Required Information LOCKED POSITION (Must be completed)**

Set hinges to LOCK at cast position  
\_\_\_\_\_ **OR** \_\_\_\_\_

Set hinge at:  0  5  10  15  Other \_\_\_\_\_

**PLEASE FOLLOW STEP-BY-STEP  
CAST PROTOCOL INSTRUCTIONS**

**Locking Joint Options**

- Single Pivot With No Free Motion (*lowest profile*)
- Single Pivot With Free Motion (*requires Cables with Twist Release*)
- 5 Bar Trigger Locks With No Free Motion
- 5 Bar Trigger Locks With Free Motion
  - Manual Lock  Manual Free Motion
  - Install Extension Assist Bands/Posts

**Cable Release Options**

- Cables With Twist Release (*centered on anterior proximal thigh band*)
- Cables With Pull Up Lever
  - Lateral Side (*recommended*)  Centered On Top Band
- Set 5 Bar Flexion Stops at:
  - 15°  30°  45°  60°  75°  90°

**Becker Knee Joints** (*Townsend stocked items*)

- Modular Ring Lock Model 1402-B
- Automatic Angled Levered Lock Model 1017A
- Modular Ratchet Lock Model 1018A

**Becker 1017 and 1018 Lock Release Options**

*Townsend Twist and Lever Release System CANNOT be used with Becker knee joints*

- Bail Lock Integrated Strap System (*BLISS Model MX-003-BLIS for use on model 1017 and 1018*)
- Bend Levers As A Bail Rod
- Cut and bend model 1017 or 1018 as manual triggers

**Condylar pads**

No  Medial  Lateral  Both

**Color/Fabric Inlay**

- |                                    |   |  |
|------------------------------------|---|--|
| <input type="checkbox"/> Black     | <input type="checkbox"/> Royal Blue     | <input type="checkbox"/> Sheer Teal*     |
| <input type="checkbox"/> Beige     | <input type="checkbox"/> Green          | <input type="checkbox"/> Sheer Purple*   |
| <input type="checkbox"/> Gray      | <input type="checkbox"/> Burgundy       | <input type="checkbox"/> US Flag Fabric* |
| <input type="checkbox"/> Red       | <input type="checkbox"/> Clear Graphite | <input type="checkbox"/> Fabric          |
| <input type="checkbox"/> Navy Blue | <input type="checkbox"/> Sheer Red*     | <small>-1 yard from patient*</small>     |

**Anti-Migration Silicon Infused Strap Pads\***

**Undersleeves\***

- 18" Cotton
- 18" Neoprene
- 22" Neoprene

**Thigh Sleeves\***

- 1/16 Comfort Thigh Sleeve

\*Indicates additional charges apply