

**Account Contact Information**

Name: \_\_\_\_\_ Email: \_\_\_\_\_ Phone: \_\_\_\_\_

**Billing and Shipping**

PO# \_\_\_\_\_ Billing Account #: \_\_\_\_\_ Shipping Account #: \_\_\_\_\_

<b>Shipping Preference</b>	Billing Address: _____	Shipping Address: _____
<input type="checkbox"/> Ground	_____	_____
<input type="checkbox"/> Next Day A.M.	City: _____	City: _____
<input type="checkbox"/> Next Day P.M.	State _____	State: _____
<input type="checkbox"/> 2-Day A.M.	Zip _____	Zip: _____
<input type="checkbox"/> 2-Day P.M.		

*(If no preference is indicated, this order will be shipped 2 Day P.M.) Note: We do not ship products directly to patients.*

**Patient Information**

**Fit Date:** \_\_\_\_\_

**Patient's Last Name:** \_\_\_\_\_

**Patient's First Name:** \_\_\_\_\_

Male  Female **Age** \_\_\_\_\_

**Weight** \_\_\_\_\_ (LBS) **Height** \_\_\_\_\_ (IN)

**Leg:**  Left  Right

**TM5+ Hinges (Includes Extension Stops)**

Optional Flexion Stop Kit\*

**Options**

- Spooner Patella Attachment\*
- Motocross Knee Guard\*
- Anti-Migration Silicon Infused Strap Pads\***

**Brace Cover\***

Posterior Closure  Pull On

**Sleeves\***


- 18" Cotton  22" Neoprene
- 18" Neoprene  1/16" Comfort Thigh Sleeve
- C/S Wrap (for compression and enhanced suspension)

**Size**

Small  Medium  Large  X-Large

Sizing Reference			
	6" above	M-L width	6" below
Small	15.5" to 18.5"	3.5" to 4"	12.25" to 13.75"
Medium	18.5" to 21"	4" to 4.5"	13.25" to 15"
Large	21" to 23.5"	4.5" to 5"	14.25" to 15.75"
X-Large	23.5" to 25"	5" to 5.5"	15" to 17"

*Please provide the three measurements, below, so we can determine the most appropriate standard size BOLD brace to ship. The knee width M-L is the PRIMARY measurement because the central aluminum structure is very rigid and not easily adjusted. The leg circumferences are SECONDARY reference measurements, as the semi-rigid nylon shells will generally contour to the leg when the thigh and tibia straps are tightened. For a disproportional leg, consider using a "custom assembled" Townsend Rebel or Rebel Pro brace.*



\*Indicates additional charges apply