

Account Contact Information

Name: _____ Email: _____ Phone: _____

Billing and Shipping

PO# _____ Billing Account #: _____ Shipping Account #: _____

Shipping Preference

- Ground
- Next Day A.M.
- Next Day P.M.
- 2-Day A.M.
- 2-Day P.M.

Billing Address: _____ Shipping Address: _____

City: _____ City: _____

State _____ Zip _____ State: _____ Zip: _____

(If no preference is indicated, this order will be shipped 2 Day P.M.) Note: We do not ship products directly to patients.

Received Date

Thuasne USA's shipping department use only

Patient Information

Fit Date: _____

Patient's Last Name: _____

Patient's First Name: _____

Male Female **Age** _____

Weight _____ (LBS) **Height** _____ (IN)

Leg: Left Right

Patient's Diagnosis: _____

Comments: _____

- Anterior Band Model**
- Three Bands *(One thigh, two tibia)* Four Bands *(Two thigh, two tibia)*
- Three Bands *(Two thigh, one tibia)*

- Polio Band Model** *(Anterior thigh, posterior calf)*
- Three Bands *(One anterior thigh, two posterior calf)*
- Three Bands *(Two anterior thigh, one posterior calf)*
- Four Bands *(One proximal anterior thigh, one distal posterior thigh, one proximal calf, one distal calf band)*

Thigh Band Length **Tibia Band Length**

7 Inch 8 Inch 7 Inch 8 Inch

9 Inch Other ____" 9 Inch Other ____"

- Graphite Extension to Floor With Heel Cup**
(Must complete graphite extension to floor form)
- C/S Package*** *Single thigh band only (for dynamic compression and enhanced suspension)*
- Set hinges to LOCK at cast position

_____ **OR** _____

Set hinge at: 0 5 10 15 -5 No Medial

-10 -15 Other ____ Lateral Both

Locking Joint Options

- Single Pivot With No Free Motion *(lowest profile)*
- Single Pivot With Free Motion *(requires Cables with Twist Release)*
- 5 Bar Trigger Lock With No Free Motion
- 5 Bar Trigger Lock With Free Motion
 - Install Extension Assist Bands/Posts

Cable Release Options

- Cables with twist release Cables with pull up lever
- Mount release lateral side Mount release center band
- Set 5 Bar Flexion Stops at:
 - 15° 30° 45° 60° 75° 90°

Becker Knee Joints *(Townsend stocked items)*

- Modular Ring Lock Model 1402-B
- Automatic Angled Levered Lock Model 1017A
- Modular Ratchet Lock Model 1018A

Becker 1017 and 1018 Lock Release Options

Townsend Twist and Lever Release System CANNOT be used with Becker knee joints

- Bail Lock Integrated Strap Bend Levers As A Bail Rod System *(BLISS) Model MX-003-BLIS (for use on model 1017 and 1018)*
- Cut and bend model 1017 or 1018 as manual triggers

Finish and Color

Textured Powdercoat Finish

- Black Royal Blue
- Antique Pewter Burgundy

High Gloss Paint Finish

- Black Emerald Green Burgundy
- Royal Blue Steel Blue White
- Burnt Orange Quicksilver Beige
- Dark Violet Indy Yellow

Custom Paint Finish* - Indicate Custom Paint # _____

Anti-Migration Silicon Infused Strap Pads*

Sleeves*

- 18" Cotton 22" Neoprene
- 18" Neoprene 1/16" Comfort Thigh Sleeve
- C/S Wrap *(for compression and enhanced suspension)*

*Indicates additional charges apply