



THUASNE

ACTIVE RELIEVER
OSTEOARTHRITIS KNEE BRACE



TOWNSEND
THUASNE USA



QUINN MEDICAL
THUASNE USA

ACTIVE RELIEVER

RELIEVER OA KNEE BRACES

Mild/moderate unicompartmental osteoarthritis, articular cartilage defect repair, meniscal cartilage repair or osteochondral defect repair.

FEATURES

- Breathable spacer fabric sleeve with posterior expandable Lycra panel
- Anatomically correct roll-and-glide Townsend Motion Hinge with adjustable ROM stops
- Semi-rigid, self-contouring nylon superstructure sewn into the sleeve
- Bendable bar for customization of the fit and counterforce
- Corrective Force Strap applies maximum counterforce at mid-stance phase of gait
- Patient friendly quick release buckles
- Open anterior thigh section for easy donning and doffing

OPTIONS

- Medial or lateral OA
- Left or right hinge position

ORDERING

Size	6" Above	6" Below
XS	13 – 15 inch	10 – 12 inch
SM	15 – 18 inch	12 – 14 inch
MD	18 – 21 inch	14 – 16 inch
LG	21 – 23 inch	16 – 18 inch
XL	23 – 26 inch	18 – 20 inch
XXL	26 – 32 inch	20 – 22 inch
XXXL	29 – 32 inch	22 – 24 inch



Adaptable patella buttress pad



Interior thigh panel can be removed and trimmed for custom fitting



Ergonomically designed quick release buckles make it easier for the patient to apply and remove the brace