



THUASNE



DORSO RIGID® 35

PDAC Approved L0472

Control thoracolumbar flexion with the Dorso Rigid 35 premium hyperextension brace.

Dorso Rigid 35

Premium Hyperextension Orthosis



ARTICULATING PELVIC BAND AND STERNAL PAD

- Spring assisted pelvic band articulates or locks at 0° or 18°
- Sternal pad self-adjusts for comfortable control


EASY TO USE LEVER LOCK

- Effortless application and tightening for enhanced compliance
- Safe and secure latch closure
- Quick connect strap button with hand loop for easy application



Easily Optimize Fit And Control

- Adjustable telescoping aluminum alloy frame
- Thoracic band widens without removing screws, and adjusts angle in 7.5° increments
- Malleable frame for optimized fit
- Back pad positions horizontally or vertically

Size	Pelvis Circumference	Pubic Symphysis To Sternum Height	
SM	23 1/2 – 29 1/2 inch	12 1/2 inch	
MS	29 1/2 – 35 1/2 inch	12 1/2 inch	
MD	29 1/2 – 35 1/2 inch	13 3/4 inch	
LS	35 1/2 – 41 1/4 inch	13 3/4 inch	
LG	35 1/2 – 41 1/4 inch	15 inch	
XL	41 1/4 – 45 1/4 inch	16 1/4 inch	



Water Resistant Materials

- Brace can be immersed in water for showering or pool rehabilitation
- Stainless steel screws, aluminum frame, and plastic hardware
- Water repellent pads