



**THUASNE**

---

**ROM-R**







**THUASNE**

## ROM-R

**RECOMMENDED L-CODE L1832/L1833**


**Recommended for:** *Patients who require post-operative fixation and/or range of motion control for post-injury stabilization.*

### Features & Benefits

- ROM hinges with drop locks and tool-free adjustable flexion and extension settings
- Sturdy aluminum struts with extendable ends that telescope from 19" to 26" in ½" increments
- Easy to adjust
- Comfortable and light weight
- Malleable paddles with individual cool pads can be contoured to the leg to achieve an intimate fit, increase stabilization and improve brace suspension
- Universal size, no measurement required

### Options

- Breeze pad or full pad configuration now available



**Full Pads Now Available**



**Length adjustment**

