



**THUASNE**

# SafeLimb

Transtibial Residual  
Limb Protector


**Recommended L5450**

Providing comfort and security during post-surgical care, SafeLimb offers supportive protection of the residual limb in case of patient falls.

*Recommended for post-operative Below Knee Amputation (BKA) patients. Offers limb protection after transtibial amputation. Prevention of knee flexion contracture following transtibial amputation.*

## Features

- Unique design, materials and adaptability enhances fit, suspension and patient compliance
- Spherical compression foam cushion prevents trauma to the surgical wound
- Length and contour of the medial bars can be adjusted to fit patient
- Anterior opening and anterior straps for ease of application
- Proprietary "SupraPatellar Strap" prevents migration
- Lightweight and breathable material




*The SafeLimb's metal bars can be contoured without tools.*



Private  
Labeling  
Available



*The residual limb is protected by a cushioned plate at the distal end of the brace*



*Townsend's proprietary "SnapLock" allows for quick, tool-free length adjustment.*

# Full Shell BK Brace

## Specialty Knee Orthosis

### Recommended L1846


The Full Shell BK Brace is the custom solution designed to control knee instabilities for below knee amputees.

*Recommended for patients with mild to moderate ACL and MCL ligament injuries, meniscal cartilage injuries, ligament laxity of the knee, knee instabilities in the coronal plane, genu recurvatum, and quadricep weakness.*

### Features


- Can accommodate any size BK socket
- Option to add Extension Assist

### Hinge Options



- |   |                                       |
|---|---------------------------------------|
| 1. Aluminum or Stainless TM5+                       | 5. Automatic Angled Lever Lock 1017-A |
| 2. Townsend Original                                | 6. Modular Ratchet Lock 1018-A        |
| 3. Single Pivot With Manual Or Cable Release        | 7. Modular Ring Lock 1402-B           |
| 4. 5 Bar Trigger Locks With Manual Or Cable Release |                                       |

### Color Options



- |                |               |                  |
|----------------|---------------|------------------|
| 1. Gloss Black | 6. Green      | 11. Sheer Teal   |
| 2. Grey        | 7. Navy Blue  | 12. Sheer Purple |
| 3. Beige       | 8. Royal Blue | 13. Flag Fabric  |
| 4. Burgundy    | 9. Clear      |                  |
| 5. Red         | 10. Sheer Red |                  |

