

GenuStart™

Combining comfort and versatility, the GenuStart line of knee braces gives patients and providers the options they need for the treatment of ligament injuries.



THUASNE

GenuStart
L1820

GenuStart ROM
L1832/L1833

Clinical Options

Design Versatility

Patient Comfort



 **TOWNSEND**
THUASNE USA

4615 Shepard Street, Bakersfield, CA 93313 • 800.432.3466 • Fax: 800.798.2722 • ThuasneUSA.com

GenuStart

L1820

GenuStart ROM

L1832/L1833

Comfortable, versatile and easy to use, the GenuStart line of knee braces offers a variety of fit and treatment options designed for a superior patient experience. Available in two lengths, wrap or sleeve styles and two hinge options, the GenuStart is a powerful new solution featuring trusted Townsend hinge technology.

Open Wrap Style




Pull-on Sleeve Style

Hinge Options:

Townsend's patented TM5 + hinge replicates the natural roll and glide of the knee joint, improving function and limiting pistoning of the hinge.

Tool-less flexion and extension system with simple labeling for ease of use. Set flexion to 0°, 45°, 60° or 90° and extension to 10° or 20°. Six flexion/extension pins included with the brace.

Safely set, adjust and secure flexion/extension points with the easy hook-and-snap hinge cap.




Recommended for: Functional treatment of the knee before or after surgery, rehabilitation following mild to moderate sprains of the cruciate and/or collateral ligaments, knee ligament injuries or ruptures and instability of the knee joint.

Features

- Easy to fit, adjust and adapt with sleeve or wrap styling, available in short and long lengths.
- Patented Townsend ROM TM5 + or TM5 hinge options replicate the natural movement of the knee.
- Alternating strap pattern prevents rotation and improves suspension.
- Open patella buttress for increased stability and support.
- Anatomically shaped medial/lateral aluminum stays for superior support.

Ordering Methods

Order the GenuStart in any hinge and length combination to ensure individualized treatment and fit.



Size	1 Circumference measurement 6" above knee center	2 Circumference measurement at knee center	3 Circumference measurement 6" below knee center
XS	13 - 15 1/2"	11 3/4 - 13"	9 3/4 - 11 3/4"
SM	15 1/4 - 18 1/2"	13 - 14"	11 3/4 - 14"
MD	18 1/2 - 20 3/4"	14 - 15"	14 - 16 1/4"
LG	20 3/4 - 23 3/4"	15 - 17"	16 1/4 - 18"
XL	23 1/2 - 23 3/4"	17 - 19"	18 - 20"
XXL	26 1/3 - 29 1/2"	19 - 21"	20 - 22"
XXXL	29 1/2 - 32"	21 - 23"	22 - 24"

GENUSTART CONFIGURATIONS

STYLE	LONG 16"		SHORT 12"	
SLEEVE				
HINGE	TM5 +	TM5	TM5 +	TM5
WRAP				