



THUASNE

GenuStart

L1820

Comfortable, versatile and easy to use with a variety of fit and treatment options designed for a superior patient experience. GenuStart is a powerful new solution featuring trusted Townsend hinge technology.




GenuStart L1820

Recommended for: Functional treatment of the knee before or after surgery, rehabilitation following mild to moderate sprains of the cruciate and/or collateral ligaments, knee ligament injuries or ruptures and instability of the knee joint.

FEATURES:

TOWNSEND
TM5 Hinge



GenuStart 16" Sleeve

Double upper
and lower
strapping on
16" length

TM5 Hinge
and form-to-fit
uprights increase
support for
recovery and
injury prevention.




GenuStart 12" Wrap

Single upper
and lower
strapping on
12" length


- Easy to fit, adjust and adapt with sleeve or wrap styling, available in short and long lengths.
- Open patella buttress for increased stability and support.
- Patented Townsend TM5 hinge options replicate the natural movement of the knee.
- Alternating strap pattern prevents rotation and improves suspension.
- Anatomically shaped medial/lateral aluminum stays for superior support.

Ordering Methods:

Choose sleeve or wrap style, long or short length and size using the 3 measurement system.



Size	1 Circumference measurement 6" above knee center	2 Circumference measurement at knee center	3 Circumference measurement 6" below knee center
XS	13 - 15 1/2"	11 3/4 - 13"	9 3/4 - 11 3/4"
SM	15 1/4 - 18 1/2"	13 - 14"	11 3/4 - 14"
MD	18 1/2 - 20 3/4"	14 - 15"	14 - 16 1/4"
LG	20 3/4 - 23 3/4"	15 - 17"	16 1/4 - 18"
XL	23 1/2 - 23 3/4"	17 - 19"	18 - 20"
XXL	26 1/3 - 29 1/2"	19 - 21"	20 - 22"
XXXL	29 1/2 - 32"	21 - 23"	22 - 24"



3MS Sizing System