

Account Contact Information

Name: _____ Email: _____ Phone: _____

Billing and Shipping

PO# _____ Billing Account #: _____ Shipping Account #: _____

Shipping Preference

- Ground
- Next Day A.M.
- Next Day P.M.
- 2-Day A.M.
- 2-Day P.M.

Billing Address: _____ Shipping Address: _____

City: _____ City: _____

State _____ Zip _____ State: _____ Zip: _____

(If no preference is indicated, this order will be shipped 2 Day P.M.) Note: We do not ship products directly to patients.

Patient Information

Fit Date: _____

Patient's Last Name: _____

Patient's First Name: _____

Male Female Age _____

Weight _____ (LBS) **Height** _____ (IN)

Leg: Left Right

Brace Configuration

Rebel Rebel Pro Rebel Lite
(only made in 13" length)

Thigh Shell Length 7" anterior 8" anterior

Tibia Shell Length 6" anterior 7" posterior
 7" anterior
 8" anterior

XXX-Large Size Available 14" & Anterior Tibia Shell Only

Strapping Options CS Package* PCL Strap*

Color

Matte Finish

Black (Standard) Grey White
 Atlantic (Light Blue) Red

Satin Finish

Lemon Orange Lime
 Fuchsia Pacific (Dark Blue)

Accessories

- Flexion Stop Kit* (Extension stops included with all Rebels)
- Extension assist bands/posts* (Pro model only)
- Quick release buckles*
- Anti-migration silicon infused strap pads*
- Spooner patella stabilizing attachment*
- X-treme Guard (Patella Protector)*

Brace Cover* (Pull-on) S/M L/XL

Sleeves*

- 18" Cotton 22" Neoprene
- 18" Neoprene 1/16" Comfort Thigh Sleeve
- C/S Wrap (for compression and enhanced suspension)

Sizing

If your patient has proportional leg sizing (see sizing reference, below) or if you are ordering a brace for stock inventory please select from the size options.

- X-Small Small Medium Large X-Large
- XX-Large **Pro model only** XXX-Large **Pro model only**

Sizing Reference

	6" above	M-L width	6" below
X-Small	12.5" to 15.5"	3" to 3.5"	11" to 12.5"
Small	15.5" to 18.5"	3.5" to 4"	12.25" to 13.75"
Medium	18.5" to 21"	4" to 4.5"	13.25" to 15"
Large	21" to 23.5"	4.5" to 5"	14.25" to 15.75"
X-Large	23.5" to 25"	5" to 5.5"	15" to 17"
XX-Large	25" to 28"	5.5" to 6"	17" to 19"
XXX-Large*	28" to 31"	6" to 6.5"	19" to 21"

Measurement Data

Required for "custom-fitted" Assembly*


Custom-fitted (Fabricated from three leg measurements)*

If your patient has a disproportionate leg, Thuasne USA will customize the assembly of your patient's brace. Please provide leg measurements beside the illustration, below.

_____ Circumference 6 inches above mid-patella

_____ **Medial-Lateral Knee Width (not circumference) at knee center**

_____ Circumference 6 inches below mid-patella



Special Instructions: _____

*Indicates additional charges apply