

Contact Information

Clinician Fitter/Assistant/Tech Other: _____
 Name: _____
 Email: _____ Phone: _____

Ordering Clinician

CPO CO CP Other: _____
 Name: _____
 Email: _____ Phone: _____

Billing & Shipping

PO#:

Billing Account#: _____
 Shipping Account#: _____

Shipping Address: _____
 City: _____ State: _____ Zip: _____

Shipping Preference

Ground Next Day AM Next Day PM 2-Day AM 2-Day PM

(If no preference is indicated, this order will be shipped 2 Day PM) Note: We do not ship products directly to patients.

Patient Information

Fit Date: _____ Diagnosis: _____

Last Name: _____

First Name: _____

Male Female Weight _____ (LBS) Leg: Left

Age _____ Height _____ (IN) Right

Comments: _____

Anterior Band Model

Three Bands (One thigh, two tibia) Four Bands (Two thigh, two tibia)
 Three Bands (Two thigh, one tibia)

Polio Band Model (Anterior thigh, posterior calf)

Three Bands (One anterior thigh, two posterior calf)
 Three Bands (Two anterior thigh, one posterior calf)
 Four Bands (One proximal anterior thigh, one distal posterior thigh, one proximal calf, one distal calf band)

Thigh Band Length 7" 8" 9" Other _____"

Tibia Band Length 7" 8" 9" Other _____"

Graphite Extension to Floor With Heel Cup
(Must complete graphite extension to floor form)

C/S Package* Single thigh band only
(for dynamic compression and enhanced suspension)

Set hinges to LOCK at cast position

_____ OR _____

Set hinge at: 0° 5° 10° 15° 5° -10° -15°
 Other _____°

Condylar pads: No Medial Lateral Both

Locking Joint Options

Single Pivot With No Free Motion (lowest profile)
 Single Pivot With Free Motion (requires Cables with Twist Release)
 5 Bar Trigger Lock With No Free Motion
 5 Bar Trigger Lock With Free Motion
 Install Extension Assist Bands/Posts

Cable Release Options

Cables with twist release Cables with pull up lever
 Mount release lateral side Mount release center band
 Set 5 Bar Flexion Stops at:
 15° 30° 45° 60° 75° 90°

Becker Knee Joints (Townsend stocked items)

Modular Ring Lock Model 1402-B
 Automatic Angled Levered Lock Model 1017A
 Modular Ratchet Lock Model 1018A
 Adjustable Extension Automatic Spring Lever
 Knee Joint With No Free Motion 1007-B6S

Becker 1017 and 1018 Lock Release Options

Townsend Twist and Lever Release System CANNOT be used with Becker knee joints

Bail Lock Integrated Strap Bend Levers As A Bail Rod
 System (BLISS) Model MX-003-BLIS Cut and bend model 1017
(for use on model 1017 and 1018) or 1018 as manual triggers

Finish and Color

Textured Powdercoat Finish

Black Royal Blue
 Antique Pewter Burgundy

High Gloss Paint Finish

Black Emerald Green Burgundy
 Royal Blue Steel Blue White
 Burnt Orange Quicksilver Beige
 Dark Violet Indy Yellow
 Custom Paint Finish* - Indicate Custom Paint # _____
 Anti-Migration Silicon Infused Strap Pads*

Sleeves*

18" Cotton 22" Neoprene
 18" Neoprene 1/16" Comfort Thigh Sleeve
 C/S Wrap (for compression and enhanced suspension)

*Indicates additional charges apply

OF-018 REV. C

Received Date *Thuasne USA's shipping department use only*

Distributed by **Thuasne USA**

4615 Shepard Street, Bakersfield, CA 93313
 Tele: 800.432.3466 • Fax: 800.798.2722

ThuasneUSA.com

