

**Account Contact Information**

Name: \_\_\_\_\_ Email: \_\_\_\_\_ Phone: \_\_\_\_\_

**Billing and Shipping**

PO# \_\_\_\_\_ Billing Account #: \_\_\_\_\_ Shipping Account #: \_\_\_\_\_

**Shipping Preference**

- Ground
- Next Day A.M.
- Next Day P.M.
- 2-Day A.M.
- 2-Day P.M.

Billing Address: \_\_\_\_\_ Shipping Address: \_\_\_\_\_

City: \_\_\_\_\_ City: \_\_\_\_\_

State \_\_\_\_\_ Zip \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

*(If no preference is indicated, this order will be shipped 2 Day P.M.) Note: We do not ship products directly to patients.*

**Patient Information (If applicable)**

**Fit Date:** \_\_\_\_\_ **Patient Shoe Size:** \_\_\_\_\_

**Patient's Last Name:** \_\_\_\_\_

**Patient's First Name:** \_\_\_\_\_

**SpryStep® Original (weight limit = 250lbs)**

Right \_\_\_ XS \_\_\_ S \_\_\_ M \_\_\_ L \_\_\_ XL

Left \_\_\_ XS \_\_\_ S \_\_\_ M \_\_\_ L \_\_\_ XL

**SpryStep® Flex (weight limit = 250lbs)**

Right \_\_\_ XS \_\_\_ S \_\_\_ M \_\_\_ L \_\_\_ XL

Left \_\_\_ XS \_\_\_ S \_\_\_ M \_\_\_ L \_\_\_ XL

**SpryStep® Max (weight limit = 250lbs)**

Right \_\_\_ XS \_\_\_ S \_\_\_ M \_\_\_ L \_\_\_ XL

Left \_\_\_ XS \_\_\_ S \_\_\_ M \_\_\_ L \_\_\_ XL

**SpryStep® Plus (weight limit = 250lbs)**

Right \_\_\_ XS \_\_\_ S \_\_\_ M \_\_\_ L \_\_\_ XL

Left \_\_\_ XS \_\_\_ S \_\_\_ M \_\_\_ L \_\_\_ XL

**SpryStep® Pediatric (weight limit = 130lbs)**

Right \_\_\_ XS \_\_\_ S \_\_\_ M \_\_\_ L \_\_\_ XL

Left \_\_\_ XS \_\_\_ S \_\_\_ M \_\_\_ L \_\_\_ XL

**Accessories**


- Refurbishment kit (pad, strap, alligator tab)\*
- BOA Strap® (not available for Pediatric)\*
  - Right AFO, buckle inside / Left AFO, buckle outside
  - Left AFO, buckle inside / Right AFO, buckle outside


1. The forefoot section can be trimmed with scissors: Approximately 1 inch can be taken off the toe end, and ½ to ¾ inch can be trimmed from EACH side of the forefoot. The material that can be cut off the forefoot is tinted blue.

2. This indicates the shortest the footplate can be trimmed with scissors

3. This indicates the lowest the calf cuff can be trimmed with scissors

**A brace that has been trimmed cannot be returned for credit or replacement**

SpryStep® & SpryStep® Flex				
Size	Height	Foot Plate Length <sup>1</sup>	Original	Flex
XS	11 ¾"	8 ½ — 9 ½"		
SM	12 ⅝"	9 ⅛ — 10 ¼"		
MD	13 ⅜"	9 ⅝ — 10 ⅞"		
LG	14 ⅙"	10 ¼ — 11 ⅜"		
XL	15"	11 ⅜ — 12"		

SpryStep® Max			
Size	Height	Foot Plate Length <sup>1</sup>	
XS	13 ¼"	7 ¾ — 8 ½"	
SM	14 ½"	8 ¼ — 9"	
MD	15 ½"	9 — 10"	
LG	16 ½"	10 — 10 ¾"	
XL	16 ½"	10 ¾ — 11 ½"	

SpryStep® Plus			
Size	Height to Tib Tuberosity	Foot Plate Length <sup>1</sup>	Calf circumference 1" below fibula head
XS	13"	8 ½ — 9 ⅝"	11 ½ — 14 ½"
S	14 ¼"	9 ⅛ — 10 ¼"	12 ¾ — 15 ¾"
M	15 ½"	9 ⅝ — 10 ⅞"	13 ½ — 17 ½"
L	16 ¾"	10 ¼ — 11 ⅜"	14 ¾ — 18 ¾"
XL	18"	11 ⅜ — 12"	15 ½ — 20 ½"

SpryStep® Pediatric				
Size	Foot Length		Calf Height	
	Min <sup>2</sup>	Max	Min <sup>3</sup>	Max
XS	3 ¾"	5 ¾"	6 ¼"	7"
SM	4 ¼"	6 ¼"	7 ½"	8 ¼"
MD	5"	7"	8 ¾"	9 ½"
LG	5 ½"	7 ½"	9 ¾"	10 ½"
XL	6"	8"	11"	11 ¾"

*Calf cuff depth = 2 ½" (the width of the calf cuff can be trimmed up to 1 ½" on each side)*

\*Indicates additional charges apply