





Stat-A-Dyne™ ESP





Stat-A-Dyne™ ESP

(Elbow-Supination-Pronation)

Core Technologies

- Three modalities for tissue elongation (*static progressive, dynamic, & dual function*)
- Bi-directional hinges
- Circular hinge for terminal range of motion
- Light-weight and rigid aircraft aluminum maintains anatomical alignment

Feature Benefit

- Infinite ROM eliminates device bottom out in supination and pronation
- Elbow ROM from -45° hyperextension to 150° flexion
- 3-point lever arm
- Only device capable of treating extension, flexion, pronation and supination.
- Off the shelf and custom cuffs for maximum soft tissue capture
- Adjustable carrying angle

Indications

- Fracture & Dislocation
- ORIF
- Crush Injury
- Radial Head Fracture
- Joint Arthroplasty
- Tendon & Ligament Repair
- Burn
- Closed Manipulation
- Other Orthopedic and Neurological Conditions