



THUASNE

Stat-A-Dyne™
Ankle



Stat-A-Dyne™ Ankle



Recommended for *Fracture & Dislocation / ORIF / Crush Injury / CVA / TBI / Joint Arthroplasty / Tendon & Ligament Repair / Burn / Closed Manipulation / Other Orthopedic and Neurological Conditions*

Core Technologies

- Three modalities for tissue elongation (static progressive, dynamic, & dual function)
- Bi-directional hinges
- Circular hinge for terminal range of motion
- Light-weight and rigid aircraft aluminum maintains anatomical alignment

Feature Benefit

- Ankle ROM: Dorsiflexion 40°
Plantar flexion 50°
Inversion 40°
Eversion 30°
- Off the shelf and custom cuffs for maximum soft tissue capture
- Posterior radius prevents heel breakdown
- Bilateral use capability
- Anti-skid pad for stand-pivot-transfers
- Maintains hip in neutral alignment