




THUASNE



## Stat-A-Dyne™ Elbow

**Recommended for** *Fracture & Dislocation / ORIF / Crush Injury / Radial Head Fracture / Joint Arthroplasty / Tendon & Ligament Repair / Burn / Closed Manipulation / Other Orthopedic and Neurological Conditions*

### Core Technologies

- Three modalities for tissue elongation (static progressive, dynamic, or dual function)
- Bi-directional hinge for flexion and extension
- Circular hinge for terminal range of motion
- Light-weight and rigid aircraft aluminum maintains anatomical alignment

### Feature Benefit

- Elbow ROM from -45° hyperextension to 150° flexion
- Counterforce bracket provides 3-point lever arm
- Off the shelf and custom cuffs for maximum soft tissue capture
- Proximal cuff opens for terminal elbow flexion
- Forearm cuff crosses wrist joint to prevent point loading
- Adjustable L-bracket maximizes patient comfort
- Adjustable carrying angle