



THUASNE



## Stat-A-Dyne™ Pronation/Supination

*Recommended for Fracture & Dislocation / ORIF / Crush Injury / Radial Head Fracture / Joint Arthroplasty / Tendon & Ligament Repair / Burn / Closed Manipulation / Other Orthopedic and Neurological Conditions*

### Core Technologies

- Three modalities for tissue elongation (static progressive, dynamic, & dual function)
- Bi-directional hinge allows both pronation & supination in a single device
- Circular hinge for terminal end range of motion
- Light-weight and rigid aircraft aluminum maintains anatomical alignment

### Feature Benefit

- Infinite ROM eliminates bottom out
- Bivalved or segmented Kydex cuffs maximize soft tissue capture and comfort
- Radial-based forearm cuff assures anatomically correct rotation
- Light-weight and low profile

