



THUASNE

Stat-A-Dyne™
Shoulder



Stat-A-Dyne™ Shoulder



Recommended for Primary & Secondary Adhesive Capsulitis "Frozen Shoulder" / Rotator Cuff Repair / Burn / Reverse Total Shoulder / MUA / Post Cortisone Injection / Other Orthopedic and Neurological Conditions

Core Technologies

- Three modalities for tissue elongation (static progressive, dynamic, & dual function)
- Bi-directional hinges for internal/external rotation and abduction/adduction
- Circular hinge for terminal range of motion
- Light-weight and rigid aircraft aluminum maintains anatomical alignment

Feature Benefit

- ROM: Abduction 0° - 180°
 - Internal Rotation 0° - 100°
 - External Rotation 0° - 120°
- Capable of providing full External/Internal Rotation and Abduction/Adduction
- Easy application and comfortable to use promoting patient compliance
- Available in ambulatory and chair models

