



**THUASNE**


---

Stat-A-Dyne™  
**WHFO Slider**



# Stat-A-Dyne™ WHFO Slider

(Wrist Hand Finger Orthosis)



**Recommended for** Fracture & Dislocation / Sub-acute Tenolysis / Crush Injury / Tendon & Ligament Repair / Burn / Degloving / CRPS. Finger and/or Wrist Contracture / Closed Manipulation / Other Orthopedic and Neurological Conditions

## Core Technologies

- Three modalities for tissue elongation (static progressive, dynamic, & dual function)
- Bi-directional hinge allows both flexion and extension in a single device
- Circular hinge for terminal end range of motion
- Light-weight and rigid aircraft aluminum maintains anatomical alignment
- Slider technology prevents joint compression and point loading

## Feature Benefit

- ROM: Extension 0° - 90° | Flexion 0° - 90° | Full Composite Fist
- Combined motion therapy for wrist and fingers
- Long Extrinsic Extensor/Flexor Stretch and Tenodesis Stretch
- Malleable finger attachments for individualized finger ROM
- Bivalved or segmented Kydex cuffs maximize soft tissue capture and comfort
- Small, medium, and large gloves