



THUASNE

Stat-A-Dyne™
WHFO



Stat-A-Dyne™ WHFO

(Wrist Hand Finger Orthosis)



Recommended for Fracture & Dislocation / ORIF / Crush Injury / Tendon & Ligament Repair / Burn / Degloving / CRPS / Finger and/or Wrist Contracture / Closed Manipulation / Other Orthopedic and Neurological Conditions

Core Technologies

- Three modalities for tissue elongation (static progressive, dynamic, & dual function)
- Bi-directional hinge allows both flexion and extension in a single device
- Circular hinge for terminal range of motion
- Light-weight and rigid aircraft aluminum maintains anatomical alignment

Feature Benefit

- ROM: Extension 0° - 90° | Flexion 0° - 90° | Full Composite Fist
- Combined motion therapy for wrist and fingers
- Long Extrinsic Extensor/Flexor Stretch and Tenodesis Stretch
- Malleable finger attachments for individualized finger ROM
- Off the shelf or custom cuffs for maximum soft tissue capture
- Small, medium, and large gloves