



THUASNE

Stat-A-Dyne™
Wrist Slider



Stat-A-Dyne™ Wrist Slider



Recommended for Fracture & Dislocation / ORIF / Crush Injury / Tendon & Ligament Repair / Burn / Degloving / CRPS / Finger and/or Wrist Contracture / Closed Manipulation / Other Orthopedic and Neurological Conditions

Core Technologies

- Three modalities for tissue elongation (static progressive, dynamic, and dual function)
- Bi-directional hinge allows both flexion and extension in a single device
- Circular hinge for terminal end range of motion
- Light-weight and rigid aircraft aluminum maintains anatomical alignment
- Slider technology prevents joint compression and point loading

Feature Benefit

- ROM: Extension 0° - 90° | Flexion 0° - 90°
- Adjustable L-bracket maximizes patient comfort
- Bivalved Kydex cuffs maximize soft tissue capture and comfort
- Orthopedic and neurological hand platforms
- Light-weight and low profile

