




THUASNE

Stat-A-Dyne™
Wrist



Stat-A-Dyne™ Wrist



Recommended for Fracture & Dislocation / ORIF / Crush Injury / Tendon & Ligament Repair / Burn / Degloving / CRPS / Finger and/or Wrist Contracture / Closed Manipulation / Other Orthopedic and Neurological Conditions

Core Technologies

- Three modalities for tissue elongation (static progressive, dynamic, & dual function)
- Bi-directional hinge allows both flexion and extension in a single device
- Circular hinge for terminal range of motion
- Light-weight and rigid aircraft aluminum maintains anatomical alignment

Feature Benefit

- ROM: Extension 0° - 90° | Flexion 0° - 90°
- Bilateral capabilities
- Adjustable L-bracket maximizes patient comfort
- Off the shelf and custom cuffs for maximum soft tissue capture
- Orthopedic and neurological hand platforms
- Light-weight and low profile