

Hand & Wrist

Comfortable and easy to use solutions that protect, stabilize and promote healing of the hand, wrist and fingers.





Ligaflex Action PDAC APPROVED A4467

Recommended for: General weakness / postoperative conditions / post-traumatic conditions / sprained wrist / wrist tendinitis / weekend or torn TFCC ligaments / wrist instability.

Features

- Anatomically designed with Thuasne’s exclusive knit technology
- Silicone insert for comfortable protection of the ulnar condyle
- Smooth non-abrasive seams
- Adjustable pressure with removable, elastic form fitting straps
- Removable palmar splint for more active movement

	Size	1	2	3	4	5	6
 Wrist CIRCUM		5 ½ – 6"	6 – 6 ½"	6 ½ – 6 ¾"	6 ¾ – 7"	7 – 7 ½"	7 ½ – 8"



Ligaflex Pro (A)
PDAC APPROVED L3908

Recommended for: Light to moderate sprains / ligament or tendon related injury/trauma / osteo-neural injuries / wrist tendinitis / Carpal Tunnel Syndrome / Carpal Tunnel without surgery / postoperative immobilization.

Features

- Anatomical, conformable splint that supports the palm of the hand
- Lightweight and breathable micro-aerated 3D fabric
- Easy-to-use Quick Lacing System
- Universal left or right

	Size	1	2
 Wrist CIRCUM		5 – 6 ¾"	6 ¾ – 9"

Ligaflex Pro+ (B)
PDAC APPROVED L3807/L3809

Recommended for: Wrist sprains / thumb sprains / scaphoid fractures / tendinitis and tendon disorders of the wrist / thumb tendon disorders / inflammatory phases of rheumatic diseases / chronic post-traumatic or post-operative irritations / arthritis / wrist osteoarthritis.

Features

- Anatomical, conformable splint that supports the palm of the hand
- Lightweight and breathable micro-aerated 3D fabric
- Thumb spica for additional support
- Easy-to-use Quick Lacing System
- Available in left or right


	Size	1	2	3	4
 Wrist CIRCUM		5 – 6"	6 – 7"	7 – 7 ¾"	8 – 9"

Dynastab Dual (C)
RECOMMENDED L-CODE L3908

Recommended for: Mild sprains of the wrist (functional treatment) / tendinitis and tendon disorders of the wrist / inflammatory phases of rheumatic diseases.

Features

- Anatomical palmer stay provides wrist support
- Antibacterial-treated 3D knit for effective moisture-wicking
- Additional strap for extra stabilization
- Universal left or right

	Size	0	1	2	3
 Wrist CIRCUM		3 ½ – 4 ¾"	5 – 6 ¼"	6 ½ – 7 ¾"	7 ¾ – 9"

Ligaflex Rhizo (D)
PDAC APPROVED L3923/L3924

Recommended for: Light to moderate sprains / ligament and tendon related trauma / rhiartrosis / osteo-neural injuries / immobilization of the thumb.

Features

- Patented wire thumb frame for immobilization of the thumb joint
- Opening at the thumb column to prevent pressure on sensitive anatomy
- Antibacterial-treated 3D knit for effective moisture-wicking
- Particularly suitable for basal thumb arthritis
- Available in left or right

	Size	1	2
 Wrist CIRCUM		5 – 6 ¾"	6 ¾ – 9"



A



B



C



D


Ligaflex Classic Open (A)

PDAC APPROVED L3908

Recommended for: Light to moderate sprains / ligament or tendon related injury/trauma / osteo-neural injuries / wrist tendinitis / Carpal Tunnel Syndrome / Carpal Tunnel without surgery / postoperative immobilization.

Features

- Fully opens to accommodate swelling and allows for easy application and removal
- Conformable and X-ray friendly, radio-transparent stay in the palm region
- Available in left or right

Size	0	1	2
 Wrist CIRCUM	3 ½ – 4 ¾"	5 – 6 ¾"	6 ¾ – 9"

Ligaflex Manu (B)

RECOMMENDED L-CODE L3807

Recommended for: Mild sprains of the wrist (functional treatment) / following carpal tunnel interventions / tendinitis and tendon pathologies of the wrist / painful sequelae of injuries / mild sprains of the thumb / resumption of activities after Pouteau Colles distal radius fractures.

Features

- Conformable and X-ray friendly, radio-transparent stays in the palm and thumb region
- Lightweight and comfortable antibacterial-treated 3D knit for effective moisture-wicking
- 2 in 1 strapping for auto adjustment and support
- Available in left or right

Size	0	1	2	3	4
 Wrist CIRCUM	4 – 5"	5 – 6"	6 – 7"	7 – 7 ¾"	8 – 9"


Ligaflex® Finger (C)

RECOMMENDED L-CODE L3923/L3924

Recommended for: Long finger injuries such as sprains / dislocations / fractures (functional treatment or post-surgical support).

Features

- Four straps, two fixed and two removable, for optimized support and fit
- Wrist strap with anatomical shape to improve comfort during wear
- Curved stay covers and protects the finger while four straps stabilize the brace
- Rigidity for strict immobilization from the wrist up to the finger
- Adaptability to fit the patient's wrist, hand, and finger angulation
- Two stay positions: palmar or dorsal for optimum immobilization
- Machine washable

Size	1	2	3	4
 Tip of Injured Finger	< 5"	5 – 6 ¼"	6 ¼ – 7 ½"	7 ½ – 8 ¾"


Ligaflex® Boxer (D)

PDAC APPROVED L3807/L3809

Recommended for: Non-displaced 4th and 5th metacarpal fracture, postoperative 4th and 5th metacarpal support / use following strict plaster cast immobilization.

Features

- 3D inner breathable fabric for day-to-day comfortable wear
- Finger fabric separator to protect the immobilized finger
- Two shapeable rigid stays for a strict immobilization from the wrist to the finger
- Palmar and dorsal stays for optimum immobilization
- Full opening for easier application & positioning
- Convenient buckle loop tightening system
- Machine washable

Size	1	2	3	4
 Wrist CIRCUM	5 – 5 ¾"	6 – 6 ¾"	7 – 7 ¾"	8 – 9"