

Full Shell Locking Orthosis

Specialty Bracing

Contact Information

Clinician Fitter/Assistant/Tech Other: _____

Name: _____

Email: _____ Phone: _____

Ordering Clinician

CPO CO CP Other: _____

Name: _____

Email: _____ Phone: _____

Billing & Shipping

PO#: _____

Billing Account#: _____

Shipping Account#: _____

Shipping Address: _____

City: _____ State: _____ Zip: _____

Shipping Preference

Ground Next Day AM Next Day PM 2-Day AM 2-Day PM

(If no preference is indicated, this order will be shipped 2 Day PM) Note: We do not ship products directly to patients.

Patient Information

Fit Date: _____ Patient ID: _____

Age _____ Weight _____ (LBS) Height _____ (IN)

Male Female Leg: Left Right

Diagnosis: _____

Surgeries (type/date): _____

Please Follow Step-By-Step Cast Protocol Instructions

M-L measurement at knee center _____

CASTING PROTOCOL: 18-20 inch length; non-weight bearing, foot dorsiflexed; quad relaxed; landmarks indicated; cut off back of the leg.

Model

Standard (with anterior thigh and calf shells) Polio (anterior thigh shell, posterior calf shell)

Thigh Shell Length 7" 8" Other _____"

Tibia Shell Length 7" 8" 9" Other _____"

Single Strut KAFO
(Must complete graphite extension to floor form)

Hinge Options

Required Information LOCKED POSITION *(Must be completed)*

Set hinges to LOCK at cast position

OR

Set hinge at:
 0° 5° 10° 15° Other _____°

Locking Joint Accessory

Install Extension Assist Bands/Posts
 Flexion Stop Kit *(factory installed only)*
 15° 30° 45° 60° 75° 90°

Locking Hinges

Townsend Hinge U0027

5 Bar Locking Joint Single Pivot Locking

Lock options

Pull Up Lever
 Lateral side *(recommended)* Centered on Top Band

Manual Trigger
 Free motion *(5 bar locking only)* No free motion

Cable with twist release
 Positioned on anterior of shell Lateral side *(must select if posterior thigh shell)*
 Free motion No free motion

Becker Hinge U0027-B

Modular Ring Lock Model 1402-B
 Automatic Angled Levered Lock Model 1017A
 Modular Ratchet Lock Model 1018A
 Bend Levers As A Bail Rod

Lock Release Option

BLISS Bail Rod

Fabric Inlay

Black Royal Blue Sheer Teal*
 Beige Green Sheer Purple*
 Gray Burgundy US Flag Fabric*
 Red Clear Graphite Fabric
 Navy Blue Sheer Red* *1 yard from patient**

Options

Condylar pads: No Medial Lateral Both

Anti-Migration Silicon Infused Strap Pads*

Sleeves*

18" Cotton 22" Neoprene
 18" Neoprene 1/4" Comfort Thigh Sleeve
 C/S Wrap *(for compression and enhanced suspension)*

*Indicates additional charges apply

OF-003 REV. C

Received Date *Thuasne USA's shipping department use only*

Distributed by **Thuasne USA**
4615 Shepard Street, Bakersfield, CA 93313
Tele: 800.432.3466 • Fax: 800.798.2722

ThuasneUSA.com