

### Contact Information

Clinician    Fitter/Assistant/Tech    Other: \_\_\_\_\_  
 Name: \_\_\_\_\_  
 Email: \_\_\_\_\_ Phone: \_\_\_\_\_

### Ordering Clinician

CPO    CO    CP    Other: \_\_\_\_\_  
 Name: \_\_\_\_\_  
 Email: \_\_\_\_\_ Phone: \_\_\_\_\_

### Billing & Shipping

PO#:

Billing Account#: \_\_\_\_\_  
 Shipping Account#: \_\_\_\_\_

Shipping Address: \_\_\_\_\_  
 City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

### Shipping Preference

Ground    Next Day AM    Next Day PM    2-Day AM    2-Day PM

*(If no preference is indicated, this order will be shipped 2 Day PM) Note: We do not ship products directly to patients.*

### Patient Information

Fit Date: \_\_\_\_\_ Patient ID: \_\_\_\_\_  
 Age \_\_\_\_\_ Weight \_\_\_\_\_ (LBS) Height \_\_\_\_\_ (IN)  
 Male    Female   Leg:  Left    Right  
 Ligament:    ACL    PCL    LCL    MCL  
 Meniscus Damage:    Medial    Lateral  
 Surgeries (type/date): \_\_\_\_\_

**Please Follow Step-By-Step Cast Protocol Instructions**

**M-L measurement at knee center** \_\_\_\_\_  
*Casting Protocol: 18-20 inch length; non-weight bearing, foot dorsiflexed; quad relaxed; landmarks indicated; cut off back of the leg.*

### Model

- Premier Reliever U0011 (dual hinge knee brace)
- Premier Reliever1 U0012 (single hinge knee brace)

### Compartment

- Unload Medial    Unload Lateral    Dual Loadshifters

Thigh Shell Length    7"    8"

Tibia Shell Length    7"    8"

### Tibia

- ① Anterior Single Band
- ② Posterior Single Band
- ③ Double Band\* TDO10 (anterior only and not available Reliever1)
- Single Strut KAFO With Heel Cup\*  
*(Must Complete Additional Form For Brace Extension)*



### Premier Reliever

*(fabricated with TM5+ Hinges and includes an adjustable extension stop kit)*

- Flexion Stop Kit\*
- Add extension assist bands/posts\*

### Hinge Material Options (dual hinge Premier Reliever only)

- 6061 Aluminum (required for adjustable correction LOADSHIFTER)
- Stainless Steel (brace will be fabricated without LOADSHIFTER)

### Premier Reliever1

*(fabricated with an aluminum TM6 Hinge and includes an adjustable extension stop kit).*

- Optional Flexion Stop Kit

### Brace Color (Select one)

### Textured Powdercoat Finish (Lightest, Most Durable Finish)

- Black    Royal Blue
- Antique Pewter (Silver)    Burgundy

### Paint Finish

- High Gloss    Beige    Dark Violet
- Black    Emerald Green    Steel Blue
- Royal Blue    White    Indy Yellow
- Burgundy    Burnt Orange    Quicksilver

### Custom High Gloss Paint Finish\*

- Provide Custom Paint # \_\_\_\_\_

### Options (not available for Reliever1)

- Add Combined Instability (PCL Strap)\*
- C/S Package\* (for dynamic compression and enhanced suspension)
- No wraparound attachment of Synergistic Suspension Strap (recommended if patient has a prominent fibular head)
- Anti-Migration Silicon Infused Strap Pads\*
- Spooner Patella Stabilizing Attachment\*

### Brace Cover\*

- Pull On

### Sleeves\*

- 18" Cotton    22" Neoprene
- 18" Neoprene    1/4" Comfort Thigh Sleeve
- C/S Wrap (for compression and enhanced suspension)

Special Instructions: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\*Indicates additional charges apply

OF-028 REV C

**Received Date** *Thuasne USA's shipping department use only*