

Contact Information

Clinician Fitter/Assistant/Tech Other: _____

Name: _____

Email: _____ Phone: _____

Billing & Shipping

PO#:

Billing Account#: _____

Shipping Account#: _____

Ordering Clinician

CPO CO CP Other: _____

Name: _____

Email: _____ Phone: _____

Shipping Address: _____

City: _____ State: _____ Zip: _____

Shipping Preference

Ground Next Day AM Next Day PM 2-Day AM 2-Day PM

(If no preference is indicated, this order will be shipped 2 Day PM) Note: We do not ship products directly to patients.

Patient Information

Fit Date: _____ Patient ID: _____

Age _____ Weight _____ (LBS) Height _____ (IN)

Male Female Leg: Left Right

Diagnosis: _____

Surgeries (type/date): _____

Please Follow Step-By-Step Cast Protocol Instructions

M-L measurement at knee center _____

Casting Protocol: 18-20 inch length; non-weight bearing, foot dorsiflexed; quad relaxed; landmarks indicated; cut off back of the leg.

Model

Anterior Band Model



Polio Band Model



Thigh Band Length 7" 8" 9" Other _____"

Tibia Band Length 7" 8" 9" Other _____"

Single Strut KAFO with Heel Cup
(Must complete graphite extension to floor form)

C/S Package* *Only available on option 1 (for dynamic compression and enhanced suspension)*

Hinge Options

Set hinges to LOCK at cast position

_____ OR _____

Set hinge at:

0° 5° 10° 15° 5° 10° 15°

Other _____°

Townsend Hinge U0013

5 Bar Locking Joint Single Pivot Locking

Lock options

Pull Up Lever

Lateral side *(recommended)* Centered on Top Band

Manual Trigger

Free motion *(5 bar locking only)* No free motion

Cable with twist release

Positioned on anterior of shell Lateral side
(must select if posterior thigh shell)

Free motion No free motion

Becker Hinge U0013-B

- Modular Ring Lock Model 1402-B
- Automatic Angled Levered Lock Model 1017A
- Modular Ratchet Lock Model 1018A
- Bend Levers As A Bail Rod

Lock Release Option

BLISS Bail Rod

Townsend Joints Accessory

- Install Extension assist Bands/Posts
 - Extension Stop Kit *(5 Bar Free & TM5 only)*
 - Flexion Stop Kit *(TM5 only)*
 - Flexion Stop for 5 Bar Free and 5 Bar Locking *(Factory installed only)*
- 15° 30° 45° 60° 75° 90°

Brace Color *(Select one)*

Textured Powdercoat Finish *(Lightest, Most Durable Finish)*

Black Royal Blue

Antique Pewter (Silver) Burgundy

Paint Finish

High Gloss Beige Dark Violet

Black Emerald Green Steel Blue

Royal Blue White Indy Yellow

Burgundy Burnt Orange Quicksilver

Custom High Gloss Paint Finish*

Provide Custom Paint # _____

Options

Condylar pads: No Medial Lateral Both

Anti-Migration Silicon Infused Strap Pads*

Sleeves* 18" Cotton 18" Neoprene 22" Neoprene

1/6" Comfort Thigh Sleeve C/S Wrap
(for compression & enhanced suspension)

*Indicates additional charges apply

OF-018 REV. D

Received Date *Thuasne USA's shipping department use only*