




THUASNE

---

# **SLEEQ<sup>®</sup> ADAPT** 639

Comfort and support  
that's truly customizable.  
The SLEEQ Adapt is the  
premier modular LSO by  
Thuasne USA.




Step-Down System






### Maximum Customization

Gives clinicians the ability to customize rigidity for support and individual comfort.



### Patented Invis-Adjust


Universal sizing technology fits 26" to 65" waist circumference



### Dual-Zone Technology

Non-elastic compression with pressure relief at the ribs and hips, even while seated!

**Recommended** for low back support, sagittal plane motion restriction (and coronal plane motion restriction on 637 & 639), increased intracavity pressure to reduce spine and disc load, acute and/or chronic low back and lower extremity pain (sciatica), and recovery from post-operative spinal procedures or injury.



### Modifiable Support

Interior support panels can be trimmed and/or heat modified for maximum control



### SLEEQ ADAPT 627

Low-profile lumbar orthosis



### SLEEQ ADAPT 631

Lumbar-sacral orthosis



### SLEEQ ADAPT 637

Lumbar-sacral orthosis with sagittal-coronal control



### SLEEQ ADAPT 639

Lumbar-sacral orthosis with sagittal-coronal control & enhanced anterior panel



### SLEEQ ADAPT EXT

Extension panel for waist accommodates 57"– 65" circumference

Part Number	Product Name	HCPC Code
901700	SLEEQ Adapt 627	L0627/L0642
901800	SLEEQ Adapt 631	L0631/L0648
901900	SLEEQ Adapt 637	L0637/L0650
902000	SLEEQ Adapt 639	L0639/L0651
901500	SLEEQ Adapt EXT	N/A

