



THUASNE

Air Townsend & Air Lite

Ligament Bracing Solutions

Contact Information

Clinician Fitter/Assistant/Tech Other: _____

Name: _____

Email: _____ Phone: _____

Ordering Clinician

CPO CO CP Other: _____

Name: _____

Email: _____ Phone: _____

Billing & Shipping

PO#: _____

Billing Account#: _____

Shipping Account#: _____

Shipping Address: _____

City: _____ State: _____ Zip: _____

Shipping Preference

Ground

Next Day AM

Next Day PM

2-Day AM

2-Day PM

(If no preference is indicated, this order will be shipped 2 Day PM) Note: We do not ship products directly to patients.

Patient Information

Fit Date: _____ Patient ID: _____

Age _____ Weight _____ (LBS) Height _____ (IN)

Male Female Leg: Left Right

Ligament: ACL PCL LCL MCL

Meniscus Damage: Medial Lateral

Surgeries (type/date): _____

Please Follow Step-By-Step Cast Protocol Instructions

M-L measurement at knee center _____

Casting Protocol: 18-20 inch length; non-weight bearing, foot dorsiflexed; quad relaxed; landmarks indicated; cut off back of the leg.

Model AirTownsend U0020 Air Lite U0020-AL

Configuration

ACL Combined Instabilities (PCL)*
 PCL Strap* PCL Rigid Band*

Thigh Shell Length 7" 8" Other_____"

Thigh Band Width 1½" 2" 2½" Other_____"

Tibia Shell Length 7" 8" 9" Other_____"

Single Strut KAFO With Heel Cup
(Must Complete Additional Form For Brace Extension)

Hinges

TM5+ Hinges – Includes extension stop kit
 Flexion stop kit*
 Add extension assist bands/posts*

Hinge Material

6061 Aluminum Stainless Steel*
(standard, if no hinge material is indicated)

Tibia Shell Anti-Rotation Bolster (Air Townsend only)

Standard Minimal Soft No Bolster

Special Trim Lines (Air Townsend only)

- Full Figure: For additional soft tissue containment
- No tibia shell hole
- Ski boot cut: ½ inch notched in distal tibia shell
- Rodeo: Medial thigh cut-out for rider comfort
- Customized Shell Design (include instructions/drawing)

Synergistic Suspension Strap Attachment

- Standard attachment (Lateral end recessed inside tibia shell)
- Lateral end attached to outer shell (for prominent fibular head)

Fabric Inlay

- Black Royal Blue Sheer Purple*
- Beige Green US Flag Fabric*
- Gray Burgundy Fabric
1 yard from patient*
- Red Clear Graphite
- Navy Blue Sheer Red*

Options

- Anti-Migration Silicon Infused Strap Pads*
- Spooner Patella Stabilizing Attachment*
- X-treme Guard (Patella Protector)*

Brace Cover* Pull On

Sleeves*

- 18" Cotton 22" Neoprene
- 18" Neoprene ¼" Comfort Thigh Sleeve
- C/S Wrap (for compression and enhanced suspension)

Special Instructions: _____

*Indicates additional charges apply

OF-002REV E

Received Date *Thuasne USA's shipping department use only*

Distributed by **Thuasne USA**

4615 Shepard Street, Bakersfield, CA 93313
Tele: 800.432.3466 • Fax: 800.798.2722

ThuasneUSA.com