



THUASNE

# Full Shell Locking Orthosis

Specialty Bracing

## Contact Information

Clinician  Fitter/Assistant/Tech  Other: \_\_\_\_\_

Name: \_\_\_\_\_

Email: \_\_\_\_\_ Phone: \_\_\_\_\_

## Ordering Clinician

CPO  CO  CP  Other: \_\_\_\_\_

Name: \_\_\_\_\_

Email: \_\_\_\_\_ Phone: \_\_\_\_\_

## Billing & Shipping

PO#:

Billing Account#: \_\_\_\_\_

Shipping Account#: \_\_\_\_\_

Shipping Address: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

## Shipping Preference

Ground

Next Day AM

Next Day PM

2-Day AM

2-Day PM

(If no preference is indicated, this order will be shipped 2 Day PM) Note: We do not ship products directly to patients.

## Patient Information

Fit Date: \_\_\_\_\_ Patient ID: \_\_\_\_\_

Age \_\_\_\_\_ Weight \_\_\_\_\_ (LBS) Height \_\_\_\_\_ (IN)

Male  Female

Leg:  Left  Right

Diagnosis: \_\_\_\_\_

Surgeries (type/date): \_\_\_\_\_

**Please Follow Step-By-Step Cast Protocol Instructions**

## M-L measurement at knee center \_\_\_\_\_

Casting Protocol: 18-20 inch length: non-weight bearing, foot dorsiflexed; quad relaxed; landmarks indicated; cut off back of the leg.

## Model

Standard (with anterior thigh and calf shells)

Polio (anterior thigh shell, posterior calf shell)

Thigh Shell Length  7"  8"  Other: \_\_\_\_\_"

Tibia Shell Length  7"  8"  9"  Other: \_\_\_\_\_"

Single Strut KAFO

(Must complete graphite extension to floor form)

## Hinge Options

**Required Information LOCKED POSITION (Must be completed)**

Set hinges to LOCK at cast position

OR

Set hinge at:

0°  5°  10°  15°  Other: \_\_\_\_\_°

## Locking Joint Accessory

Install Extension Assist Bands/Posts

Flexion Stop Kit (factory installed only)

15°  30°  45°  60°  75°  90°

## Locking Hinges

### Townsend Hinge U0027

5 Bar Locking Joint

Single Pivot Locking

## Lock Release Option

Manual Trigger

Cable w/Twist release on anterior side

Cable w/Twist release on lateral side (if posterior configuration - has to be lateral)

## Locking Motion

Free Motion

(not available with lever)

No Free Motion

## Becker Hinge U0027-B

Modular Ring Lock Model 1402-B

Automatic Angled Levered Lock Model 1017A

Modular Ratchet Lock Model 1018A

Bend Levers As A Bail Rod

## Lock Release Option

BLISS

Bail Rod

## Fabric Inlay

Black

Beige

Gray

Red

Navy Blue

Royal Blue

Green

Burgundy

Clear Graphite

Sheer Red\*

Sheer Purple\*

US Flag Fabric\*

Fabric

1 yard from patient\*

## Options

**Condylar pads:**  No

Medial

Lateral

Both

Anti-Migration Silicon Infused Strap Pads\*

## Sleeves\*

18" Cotton

18" Neoprene

C/S Wrap (for compression and enhanced suspension)

22" Neoprene

1/8" Comfort Thigh Sleeve

\*Indicates additional charges apply

OF-003 REV D

Received Date *Thuasne USA's shipping department use only*

Distributed by **Thuasne USA**

4615 Shepard Street, Bakersfield, CA 93313

Tele: 800.432.3466 • Fax: 800.798.2722

ThuasneUSA.com