



THUASNE

# Full Shell Series

## Osteoarthritis Bracing Solutions

### Contact Information

Clinician  Fitter/Assistant/Tech  Other: \_\_\_\_\_

Name: \_\_\_\_\_

Email: \_\_\_\_\_ Phone: \_\_\_\_\_

### Ordering Clinician

CPO  CO  CP  Other: \_\_\_\_\_

Name: \_\_\_\_\_

Email: \_\_\_\_\_ Phone: \_\_\_\_\_

### Billing & Shipping

PO#: \_\_\_\_\_

Billing Account#: \_\_\_\_\_

Shipping Account#: \_\_\_\_\_

Shipping Address: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

### Shipping Preference

Ground

Next Day AM

Next Day PM

2-Day AM

2-Day PM

*(If no preference is indicated, this order will be shipped 2 Day PM) Note: We do not ship products directly to patients.*

### Patient Information

Fit Date: \_\_\_\_\_ Patient ID: \_\_\_\_\_

Age \_\_\_\_\_ Weight \_\_\_\_\_ (LBS) Height \_\_\_\_\_ (IN)

Male  Female Leg:  Left  Right

Ligament:  ACL  PCL  LCL  MCL

Meniscus Damage:  Medial  Lateral

Surgeries (type/date): \_\_\_\_\_

**Please Follow Step-By-Step Cast Protocol Instructions**

### M-L measurement at knee center \_\_\_\_\_

*Casting Protocol: 18-20 inch length: non-weight bearing, foot dorsiflexed; quad relaxed; landmarks indicated; cut off back of the leg.*

### Model

Reliever U0025  Reliever Air U0025-RA  Reliever Air Lite U0025-RAL

### Compartment

Unload Medial  Dual Loadshifters  Unload Lateral

Thigh Shell Length  7"  8"  9"

Tibia Shell Length  7"  8"  9"  Other \_\_\_\_\_"

Single Strut KAFO With Heel Cup\*  
*(Must Complete Additional Form For Brace Extension)*

**Hinges** *Brace is fabricated with LOADSHIFTER (no charge) The LOADSHIFTER can be used to increase correction by shifting the thigh shell angle. ONLY offered with TM5+ Aluminum Hinges.*

### Hinge Material

6061 Aluminum *(required for LOADSHIFTER)*  
 Stainless Steel\* *(brace will be fabricated without LOADSHIFTER)*

### Accessory

Flexion stop kit\*  
 Add extension assist bands/posts\*

### Tibia Shell Anti-Rotation Bolster *(Air Townsend only)*

Standard  Minimal  Soft  No Bolster

### Special Trim Lines *(Air Townsend only)*

Full Figure: *(Reliever Air Only):* Helps contain soft tissue  
 Customized Shell Design *(include instructions/drawing)*

### Synergistic Suspension Strap Attachment

Standard attachment *(Lateral end recessed inside tibia shell)*  
 Lateral end attached to outer shell *(for prominent fibular head)*  
 Double Rivet Suspension Strap

### Fabric Inlay

Black  Royal Blue  Sheer Purple\*  
 Beige  Green  US Flag Fabric\*  
 Gray  Burgundy  Fabric  
 Red  Clear Graphite *1 yard from patient\**  
 Navy Blue  Sheer Red\*

### Options

Anti-Migration Silicon Infused Strap Pads\*

Spooner Patella Stabilizing Attachment\*

### Brace Cover\*

Pull On

### Sleeves\*

18" Cotton  22" Neoprene  
 18" Neoprene  1/4" Comfort Thigh Sleeve  
 C/S Wrap *(for compression and enhanced suspension)*

Special Instructions: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\*Indicates additional charges apply

OF-023 REV D

Received Date *Thuasne USA's shipping department use only*

Distributed by **Thuasne USA**

4615 Shepard Street, Bakersfield, CA 93313  
Tele: 800.432.3466 • Fax: 800.798.2722

ThuasneUSA.com