

**Account Contact Information**

Name: \_\_\_\_\_ Email: \_\_\_\_\_ Phone: \_\_\_\_\_

**Billing and Shipping**

PO# \_\_\_\_\_ Billing Account #: \_\_\_\_\_ Shipping Account #: \_\_\_\_\_

Billing Address: \_\_\_\_\_ Shipping Address: \_\_\_\_\_

- Shipping Preference**
- Ground
  - Next Day A.M.
  - Next Day P.M.
  - 2-Day A.M.
  - 2-Day P.M.

City: \_\_\_\_\_ City: \_\_\_\_\_

State \_\_\_\_\_ Zip \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

*(If no preference is indicated, this order will be shipped 2 Day P.M.) Note: We do not ship products directly to patients.*

**Patient Information**

**Patient's Last Name:** \_\_\_\_\_

**Patient's First Name:** \_\_\_\_\_

**Measuring Information**

Measure arm from elbow crease to base of the index finger:  
S/M up to 14 inches; L/XL up to 15 ½ inches.

**Low Profile ROM Elbow (Dual Hinges) With Shells And Optional Hand Extension**

- Flexion can be limited at 90, 110 and 120 degrees, and extension can be stopped at 0, 15, 30, 45 and 90 degrees.
- The screw stops are adjusted using a small flathead screwdriver.
- The hinge can be locked at 90 degrees.

Quantity: \_\_\_\_\_ S/M \_\_\_\_\_ L/XL

Install Optional Hand Extension? \_\_\_\_\_ Yes \_\_\_\_\_ No

Include Shoulder Sling Strap? \* \_\_\_\_\_ Yes \_\_\_\_\_ No

**Multi-Adjustment ROM Elbow**

- Comes with four straps and individual black pads that fit under each strap.
- Tool-less flexion and extension stops are in 10 degree increments.
- Locking is possible by offsetting the flexion and extension stops at the same degree.

**Single Hinge (SH) Model**

Quantity: \_\_\_\_\_ LEFT ARM Length: \_\_\_\_\_ 15 inch \_\_\_\_\_ 16 inch

Quantity: \_\_\_\_\_ RIGHT ARM Length: \_\_\_\_\_ 15 inch \_\_\_\_\_ 16 inch

**Dual Hinge (DH) Model**

Length (Indicate quantity): \_\_\_\_\_ 15 inch \_\_\_\_\_ 16 inch

**ROM-LBO**

- The LBO comes with dual or single uprights, four straps, tool less adjustability, and 8 individual Breeze pads.
- Tool-less flexion and extension stops are in 10 degree increments.
- The LBO can be locked at 90 degrees for post-surgical protocols.
- Adjustable length 14 ½" – 17"

**Single Hinge (SH) Model**

Quantity: \_\_\_\_\_ LEFT ARM

Quantity: \_\_\_\_\_ RIGHT ARM

**Dual Hinge (DH) Model**

Quantity: \_\_\_\_\_ Universal

\*Indicates additional charges apply