



**THUASNE®**

# Select Walker Series

PDAC APPROVED L4360/L4361

NON-PNEUMATIC L4386/L4387

**Recommended for:** *Grade 1, 2 and 3 ankle sprains and stable or stable fixed fractures of the foot or ankle.*

## Features

- Offered in
  - > **Fixed ankle**
  - > **Fixed ankle with air bootie**
- Formable plastic uprights
- Comfortable padded insole
- Shock absorbing outer sole rocker allows natural gait
- Foam bootie available with or without air bladder

## Options

- Full height or mid-calf uprights
- Heel wedges (optional)
- Universal left or right

## Sizing

Size	XS	SM	MD	LG	XL
<b>Men US Shoe Size</b>	2-4	4-7	7-10	10-13	13-16
<b>Women US Shoe Size Height</b>	3½-5	5-8	8-11	11-15	15+



**Mid-Calf Uprights**



**Full Height Uprights**



### Anatomically Shaped Uprights

*Provides rigid support for security and stability.*



### Integrated Strapping

*Helps aid in full circumferential support and compression.*



### Shock Absorbing Sole

*Low profile allows for normal gait and reduces the force upon heel strike.*