

Contact Information

Clinician Fitter/Assistant/Tech Other: _____

Name: _____

Email: _____ Phone: _____

Billing & Shipping

PO#: _____

Billing Account#: _____

Shipping Account#: _____

Ordering Clinician

CPO CO CP Other: _____

Name: _____

Email: _____ Phone: _____

Shipping Address: _____

City: _____ State: _____ Zip: _____

Shipping Preference

Ground Next Day AM Next Day PM 2-Day AM 2-Day PM

(If no preference is indicated, this order will be shipped 2 Day PM) Note: We do not ship products directly to patients.

Patient Information (If applicable)

Fit Date: _____ Patient ID: _____

Shoe Size: _____

SpryStep® Original (weight limit = 250lbs)

Right ____ XS ____ S ____ M ____ L ____ XL

Left ____ XS ____ S ____ M ____ L ____ XL

SpryStep® Flex (weight limit = 250lbs)

Right ____ XS ____ S ____ M ____ L ____ XL

Left ____ XS ____ S ____ M ____ L ____ XL

SpryStep® Max (weight limit = 250lbs)

Right ____ XS ____ S ____ M ____ L ____ XL

Left ____ XS ____ S ____ M ____ L ____ XL

SpryStep® Plus (weight limit = 250lbs)

Right ____ XS ____ S ____ M ____ L ____ XL

Left ____ XS ____ S ____ M ____ L ____ XL

SpryStep® Pediatric (weight limit = 130lbs)

Right ____ XS ____ S ____ M ____ L ____ XL

Left ____ XS ____ S ____ M ____ L ____ XL

Accessories

- Refurbishment kit (pad, strap, alligator tab)*
- BOA Strap® (not available for Pediatric)*
 - Right AFO, buckle inside / Left AFO, buckle outside
 - Left AFO, buckle inside / Right AFO, buckle outside

1. The forefoot section can be trimmed with scissors: Approximately 1 inch can be taken off the toe end, and ½ to ¼ inch can be trimmed from EACH side of the forefoot. The material that can be cut off the forefoot is tinted blue.
2. This indicates the shortest the footplate can be trimmed with scissors
3. This indicates the lowest the calf cuff can be trimmed with scissors

A brace that has been trimmed cannot be returned for credit or replacement

Original Flex Size Height Foot Plate Length 1

		XS	11 ¼"	8 ½ – 9 ⅝"
		SM	12 ⅝"	9 ⅙ – 10 ¼"
		MD	13 ⅜"	9 ⅝ – 10 ⅞"
		LG	14 ⅙"	10 ¼ – 11 ⅜"
		XL	15"	11 ⅜ – 12"

Max Size Height Foot Plate Length 1

	XS	13 ¼"	7 ¾ – 8 ½"
	SM	14 ½"	8 ¼ – 9"
	MD	15 ½"	9 – 10"
	LG	16 ½"	10 – 10 ¾"
	XL	16 ½"	10 ¾ – 11 ½"

Plus Size Height to Tib Tuberosity Foot Plate Length 1 Calf circumference 1" below fibula head

	XS	13"	8 ½ – 9 ⅝"	11 ½ – 14 ½"
	S	14 ¼"	9 ⅙ – 10 ¼"	12 ¾ – 15 ¾"
	M	15 ½"	9 ⅝ – 10 ⅞"	13 ½ – 17 ½"
	L	16 ¾"	10 ¼ – 11 ⅜"	14 ¾ – 18 ¾"
	XL	18"	11 ⅜ – 12"	15 ½ – 20 ½"

Pediatric Size Foot Length Calf Height

		Min ²	Max	Min ³	Max
	XS	3 ¾"	5 ¾"	6 ¼"	7"
	SM	4 ¼"	6 ¼"	7 ½"	8 ¼"
	MD	5"	7"	8 ¾"	9 ½"
	LG	5 ½"	7 ½"	9 ¾"	10 ½"
	XL	6"	8"	11"	11 ¾"

Calf cuff depth = 2 ½" (the width of the calf cuff can be trimmed up to 1 ½" on each side)

*Indicates additional charges apply

OF-030 REV D

Received Date Thuasne USA's shipping department use only