



THUASNE

Select Reliever

PDAC Approved L1843/ L1851



Single upright medial/lateral
Osteoarthritis Knee Brace

Select Reliever

PDAC Approved L1843/ L1851

Give your patients the osteoarthritis relief they need with the Select Reliever from Thuasne USA. Improve patient compliance and manage inventory with this low profile, universal design perfect for medial and lateral OA.



Quick release
buckle system

Features

- Lightweight, low-profile, single upright OA brace.
- Quick release buckle system for easy-on and easy-off eliminates need for strap adjustments with every use.
- Maximized 3 points of leverage for medial compartment.
- Specially designed to correct medial and lateral OA with a lateral frame application.
- Lateral uprights avoid opposite limb during gate, increasing comfort and patient compliance.




Available Accessories

Part Number	Accessories
UU031211000	OA Knee Brace Allen Key
UU031212000	OA Knee Brace Strap+Padding
UU031213001	OA Liner Full Set Right
UU031213002	OA Liner Full Set Left
UU031214000	OA Buckles
UU031215000	One Extra Long 38" Strap

Ordering

Part Number	Hinge Side
U0312129900352	Left Lateral/Right Medial
U0312129900351	Right Lateral/Left Medial



Adjustable force correction

(Images seen display right leg offloading)