

# MOBIDERM® Standard



☐ Order    ☐ Quotation    ☐ Renewal    **Date:** \_\_\_\_ / \_\_\_\_ / \_\_\_\_

**Patient Name:** \_\_\_\_\_

**Medical Professional Name:** \_\_\_\_\_


**Fitter Name:** \_\_\_\_\_

**Case No. for renewal:** \_\_\_\_\_

☐ 1st Treatment

**Phone:** \_\_\_\_\_ **Email:** \_\_\_\_\_

## Mitten & Glove



**cA:** Hand circumference  
(head of metacarpal bones)  
**cC:** Wrist circumference

☐ Left    ☐ Right

*\*Please fill one out for each side*

**cA** \_\_\_\_\_ (cm)

**cC** \_\_\_\_\_ (cm)

*\*Note: Sizing chart is the same for both the Mitten and Glove styles.*


**Comments** \_\_\_\_\_

| Size | cC<br>(in cm) | cA<br>(in cm) | <input type="checkbox"/> Mitten |                 | <input type="checkbox"/> Glove |                 |
|------|---------------|---------------|---------------------------------|-----------------|--------------------------------|-----------------|
|      |               |               | Reference                       |                 | Reference                      |                 |
|      |               |               | right/left                      |                 | right                          | left            |
| 1    | 14-16         | 15-17         | <input type="checkbox"/>        | 373101 20100299 | <input type="checkbox"/>       | 373201 20100201 |
| 2    | 16-18         | 17-19         | <input type="checkbox"/>        | 373101 20200299 | <input type="checkbox"/>       | 373201 20200201 |
| 3    | 17-19         | 19-21         | <input type="checkbox"/>        | 373101 20300299 | <input type="checkbox"/>       | 373201 20300201 |
| 4    | 19-21         | 21-23         | <input type="checkbox"/>        | 373101 20400299 | <input type="checkbox"/>       | 373201 20400201 |
| 5    | 20-22         | 23-25         | <input type="checkbox"/>        | 373101 20500299 | <input type="checkbox"/>       | 373201 20500201 |
| 6    | 21-23         | 25-27         | <input type="checkbox"/>        | 373101 20600299 | <input type="checkbox"/>       | 373201 20600201 |

**Suggested HCPCS:** A6520    **Qty** \_\_\_\_\_

**Description:** Gradient compression garment, glove, padded, for nighttime use, each

## Armsleeve



**cG:** Arm circumference  
**cE:** Elbow circumference  
**cC:** Wrist circumference  
**l:** Length from wrist to axilla

☐ Left    ☐ Right

*\*Please fill one out for each side*

**cG** \_\_\_\_\_ (cm)

**cE** \_\_\_\_\_ (cm)

**cC** \_\_\_\_\_ (cm)

**l** \_\_\_\_\_ (cm)

**Comments** \_\_\_\_\_

| Size | cC<br>(in cm) | cE<br>(in cm) | cG<br>(in cm) | Reference (right/left)   |                        |
|------|---------------|---------------|---------------|--------------------------|------------------------|
|      |               |               |               | Normal<br>l = 40 - 45 cm | Long<br>l = 45 - 50 cm |
| 1    | 14-16         | 22-25         | 25-29         | <input type="checkbox"/> | 373002 20100232        |
| 1+   |               | 24-27         | 28-32         | <input type="checkbox"/> | 373002 40100232        |
| 2    | 16-18         | 25-28         | 29-33         | <input type="checkbox"/> | 373002 20200232        |
| 2+   |               | 27-30         | 31-35         | <input type="checkbox"/> | 373002 40200232        |
| 3    | 18-20         | 28-31         | 32-36         | <input type="checkbox"/> | 373002 20300232        |
| 3+   |               | 30-33         | 35-39         | <input type="checkbox"/> | 373002 40300232        |
| 4    | 20-22         | 31-34         | 35-39         | <input type="checkbox"/> | 373002 20400232        |
| 4+   |               | 33-36         | 38-42         | <input type="checkbox"/> | 373002 40400232        |
| 5    | 22-24         | 34-37         | 38-42         | <input type="checkbox"/> | 373002 20500232        |
| 5+   |               | 36-39         | 41-45         | <input type="checkbox"/> | 373002 40500232        |

**Suggested HCPCS:** A6522    **Qty** \_\_\_\_\_

**Description:** Gradient compression garment, arm, padded, for nighttime use, each

### Fitting Tips:

- ✓ For higher upper arm coverage, it is recommended to size up on length.
- ✓ For the Mitten style, it is recommended to size up for improved MCP and thumb comfort.